



2023 **Product
Catalogue**

Homecraft® was founded in the 1940s to honour soldiers and armed forces veterans, offering dignity to those who need extra support on life's journey. Although times have changed, the fundamental values of Homecraft have remained steadfast.

Today, Homecraft fosters independence for individuals of all ages with an extensive selection of daily living aids available globally, including bathing, toileting, kitchen aids and more. Whether you are a clinician, retailer, consumer or family, Homecraft provides peace of mind to live life to the fullest.

Be your ideal self

Homecraft is at the heart of every home when you need us. Our mission is to empower individuals in their environments. Whether it is buttoning a shirt, taking a bath or using the toilet, our aids to daily living (ADLs) are there to support you.

We know daily tasks can be challenging. Homecraft was built with the spirit and belief that each carefully crafted product will empower you to lead your own life. People turn a house into a home, and Homecraft embraces your safe space with designs to fit your lifestyle.



Table of contents

When selecting a Homecraft product, you can expect to experience the difference. Homecraft represents:

Caring spirit

The Homecraft story begins with supporting war veterans, and that heart and compassion extend to all today. Seventy-five years ago, Homecraft identified a need to craft tools for one-of-a-kind situations to keep individuals self-sufficient in their homes. Homecraft is truly the heart of the home.

Dignity

Our intentionally designed products offer Homecraft customers pride and trust to navigate their daily lives. We aim to preserve individuals' confidence and dignity so they may live life fully.

Independence

Homecraft helps individuals declare and craft their own choices and freedom. With an extensive breadth of products, Homecraft has a tool for each situation, giving people the ability to maintain a grasp on life.

Peace of mind

Our attention to detail offers the right solutions to address every need, giving you and your loved ones peace of mind.

Bathing & Toileting

Bathing

Bath Lifts	2
Bath Boards and Seats	3
Bath Mats	4
Wall-mounted Shower Seats	5
Shower Stools	6
Wheeled Shower Commode Chairs	11
Bath Accessories	14
Bath Grab Bars and Rails	16

Toileting

Raised Toilet Seats	17
Support Rails	19
Toilet Frames	20
Toileting Accessories	23

Aids to Daily Living

Dressing Aids

Clothing Aids	21
Shoehorns	25
Socks and Stocking Aids	26

General Household

Home Accessories	27
------------------	----

Reachers

Pickup Reacher	29
Handi-Reacher™	30

Bed & Seating Equipment

Tables and Trays	31
Bed and Chair Raisers	33
Bed Grab and Side Rails	34
Bed Assists	35
Chairs and Stools	37

Kitchen & Dining

Queens Cutlery	39
Caring Cutlery	40
Lightweight Foam Cutlery	41
Newstead Cutlery	42
Cutlery	43
Kings Cutlery	44
Workstation	47
Plates and Boards	48
Cups and Mugs	50





Bathing & Toileting



Bathmaster™ W Bath Lift

Quickly and easily overcome difficulties getting into and out of the bath with the Bathmaster Deltis Bath Lift designed and manufactured in Germany. This bath lift separates into two parts for easy lifting and carrying. Lithium-ion batteries offer improved performance and battery life. The fully waterproof motor is designed for long-term protection against failure or damage. High-quality materials and construction offer the utmost safety and comfort with no servicing required. It is easy to assemble, and no technical knowledge is needed. Includes bath lift, blue covers, hand controller, one-hour recharger, side flaps and two side-flap protectors.

Key features:

- Ideal for long-term use
- Compact frame suitable for nearly all baths
- Simple to add and remove from the tub
- Easy to clean and corrosion-resistant

Specifications:

Seat size	50 x 37.5 cm (19¾ x 14¾") excluding side flaps
	50 x 69.5 cm (19¾ x 27¼") including side flaps
Backrest size	68 x 35 cm (26¾ x 13¾")
Seat height	6.9 - 45.5 cm (2¾ - 18")
Footprint	56 x 29.5 cm (22 x 11½")
Seat weight	6.8 kg (15 lbs)
Backrest weight	3.7 kg (8¼ lbs)
Maximum user weight	140 kg (308 lbs)
Average number of lift cycles	9 (based on 80 kg)
Warranty	2 years



Bathmaster
deltis

Part #	Description
PAT-091324086	Bathmaster Deltis Complete without Covers
PAT-091324169	Bathmaster Deltis Complete with Blue Covers
PAT-091324557	Bathmaster Deltis Blue Covers
PAT-091324532	Bathmaster Deltis Blue Swivel Seat

Savanah™ Slatted Bath Board

The Savanah Slatted Bath Board has a robust, comfortable, quick-draining design featuring smooth, reinforced plastic slats. Rubber pads under the slats help prevent the bath board from slipping on the bath rim. Brackets can be angled to position the board further into the bath to help secure it.

Key features:

- Stainless steel, corrosion-resistant fittings
- Brackets include a serrated face and rubber buffers
- Slats fully sealed against water ingress

Part #	Length x Width	Weight
PAT-AA1091A	66 x 23.2 cm (26" x 9¼")	2.0 kg (4¼ lbs)
PAT-AA1091B	71.1 x 23.2 cm (28" x 9¼")	2.2 kg (4⅘ lbs)
PAT-AA1091C	76.2 x 23.2 cm (30" x 9¼")	2.3 kg (5 lbs)



Savanah Slatted Shower Board

The Savanah Slatted Shower Board is similar in function and appearance to the Savanah Slatted Bath Board but supplied with two additional slats for extra comfort and support. It offers a large surface area for users who need to shower on the board and extra width to give users more confidence when transferring. Savanah Slatted Shower Boards may also be used with Savanah Slatted Bath Board Handle.

Key features:

- Extra width gives users more confidence when transferring
- Stainless steel, corrosion-resistant fittings
- Made with high-quality, durable materials

Part #	Length x Width	Weight
PAT-AA1093	68.6 x 35.6 cm (27" x 14")	3.0 kg (6¾ lbs)
PAT-AA1093B	71.1 x 35.6 cm (28" x 14")	3.2 kg (7 lbs)



Savanah Slatted Bath / Shower Board Handle

The Savanah Slatted Bath Board Handle is an optional accessory for the bath or shower board to ease transfer and offer additional user security. The moulding is covered with a comfortable, soft sleeve to improve handle grip. Fitting the handle may remove the need for a wall-mounted grab rail in many circumstances.

Part #	Description	Weight
PAT-AA1092	Savanah Slatted Bath Board Handle	0.2 kg (7 oz)



Savanah Modular Bath Step

This attractive, modular bath step has been specifically designed to make daily activities, especially getting into and out of the bath, easier. Supplied as a single step with a 100mm raise, the modularity of the Savanah Bath Step allows it to be stacked on top of each other or clipped together to give a larger area or create a small stair.

Specifications:

Step raise	10 cm (4")
Overall surface area	35.6 - 45.7 cm (14 - 18")
Footprint	38.1 - 48.3 cm (15 - 19")
Additional step	5 cm (2")
Weight	1.3 kg (2¾ lb)



Part #	Description
PAT-AA1824Y	Savanah Modular Bath Step

Everyday Bath Mat

The Bath and Shower Mat is a practical solution to bathing security. The basket weave surface reduces the risk of slipping while standing or sitting, and suction cup grips hold the mat in place.

Key features:

- Made from high-grade rubber
- Multiple drainage holes to drain water quickly
- Available in three sizes

Part #	Size	Dimensions
PAT-AA1802AY	Regular	57 x 35 cm (22½ x 13¾")
PAT-AA1804A	Large	76 x 35 cm (30 x 13¾")
PAT-AA1805	Extra Large	94.5 x 35 cm (37 x 13¾")



Sure Tread Shower Mat

The non-slip Sure Tread Bath and Shower Mat gives the user total peace of mind. With raised ridges and a unique design, the mat allows the user to stand or sit without slipping on a wet surface. Enjoy the ultimate secure experience of bathing with Sure Tread.

Key features:

- Made from high-grade medical rubber
- Machine-washable at 40°C (104°F)

Part #	Description	Dimensions
PAT-AA1794	Shower Mat	55.5 x 55.5 cm (21¾ x 21¾")



Transfer Bench with Back

The Transfer Bath Bench offers an extra-wide bath seat that sits with two legs inside and two legs outside the bath and is designed to support users with limited mobility, helping increase safety during transfers. The user sits securely on the seat and has plenty of room to manoeuvre during a side transfer. Optional suction legs adds further stability and slip resistance.

Key features:

- Height adjustable, aluminium, lightweight frame
- Removable backrest
- Seat with side handle
- Aluminium frame can be reversed for use on either side of the bath
- ARTG 192448

Specifications:

Seat Width (Between Arms)	59 cm (23¼")
Seat Depth	39 cm (15¼")
Overall Width (Footprint)	79 cm (31")
Overall Depth (Footprint)	39 cm (15¼")
Seat Height Range	45 - 57 cm (17¾ - 22½")
Armrest Height (from seat)	14 cm (5½")
Weight	5.1kg (11¼ lbs)

130 kg
20½ st 286½ lbs
Maximum user weight



Part #	Description
NOV-AC36	Transfer Bench with Back
NOV-426AA	Suction Cup Legs - Set of 2

Savanah Shower Seat

The Savanah Shower Seat is a wall-mounted shower seat that conveniently folds neatly against the wall when not used. Made from aluminium and plastic, it is corrosion-resistant and can be easily cleaned. The optional backrest includes arms that fold away independently for side transfers and can be mounted at the desired user height.

Key features:

- Generous seat size
- Durable
- Height-adjustable legs
- ARTG 210119

Specifications:

Seat width	43 cm (17")
Seat depth	36 cm (14¼")
Seat height	45.5 - 58.5 cm (18 - 23")
Distance from wall	43 cm (17")
Distance from wall when folded	13 cm (5")
Weight	3.8 kg (8 lbs)

159 kg
25 st 350½ lbs
Maximum user weight



PAT-091171511



PAT-091171511

shown with optional back and arms



Part #	Description
PAT-091171511	Savanah Shower Seat
PAT-091171529	Optional Back and Arms



Sherwood Plus Bariatric Perching and Shower Stool

The Sherwood Plus Bariatric Perching and Shower Stool is extra-wide and offers generously padded vinyl-covered foam seats for additional comfort. The cross-braced construction and 2.5 cm (1") epoxy-coated steel tubing provide ultra-strong seating support. Choose from two height-adjustable shower stools and two perching stools.

Key features:

- Height-adjustable legs with pin clips and non-slip feet
- Shower stool available with depth-adjustable backrest
- Perching stools feature a slightly lower frame at the front to give the seat a 5.1-cm (2") incline
- ARTG 210119

Specifications:

Shower stool dimensions:

Seat width	61 cm (24")
Seat depth	43 cm (17")
Height adjustment	50.8 - 68.6 cm (20 - 27")
Width between arms	67 cm (26½")

255 kg
40 st 563¼ lbs
Maximum user weight

Part #	Description
PAT-091156686	Shower Stool
PAT-091156694	Shower Stool with Adjustable Padded Back



Lightweight Kitchen Stool

This lightweight aluminium kitchen stool is height adjustable, rust resistant and so can be used in the shower, as well as your homecare environment. White padded seat and back.

Specifications:

Width between arms	45.5 cm (18")
Seat depth	30 cm (11¾")
Seat to top of backrest	30 cm (11¾")
Leg height adjustment	50 cm (19¾")
Floor to seat	63.5 cm (25")

100 kg
15¼ st 220½ lbs
Maximum user weight

Part #	Description
HOM-549B	Lightweight Kitchen Stool



White Line Folding Shower Stool and Chair

The Homecraft White Line Shower Stool and Chair has been cleverly designed for optimizing space. The stool and chair both fold down to allow more space when not in use or for transporting. Height-adjustable aluminium frames with moulded plastic seats featuring a cut-out for easier personal cleansing. The modern design frame is lightweight and corrosion-resistant. Supplied flat-packed, the stool and chair is easy to assemble by unfolding and clicking into place.

Key features:

- Lightweight, aluminum frame with moulded plastic seat
- Height-adjustable legs with non-slip feet
- Folding option helps with storage and transportation
- Side arms to make it easier getting on and off
- ARTG 174421



JAN-BE02



JAN-BE04

Specifications:

Overall height	66 - 76 cm (26 - 30")
Overall width	47 cm (18½")
Overall depth	40 cm (15¾")
Seat height adjustment	47 - 57 cm (18½ - 22½")
Seat width	40 cm (15¾")
Seat depth	40 cm (15¾")

Specifications:

Overall height	88 - 98 cm (34½ - 38½")
Overall width	44 cm (17¼")
Overall depth	42 cm (16½")
Seat height adjustment	47 - 57 cm (18½ - 22½")
Seat width	40 cm (15¾")
Seat depth	40 cm (15¾")
No. of adjustment heights	5

Part #	Description	Ext. dimensions (W x D)
JAN-BE02	Shower Stool with Arms	47 x 41 cm (18½ x 16")
JAN-BE04TP001W	Rubber Tip	

Part #	Description	Ext. dimensions (W x D)
JAN-BE04	Folding Shower Chair	44 x 42 cm (17¼ x 16½")
JAN-BE04TP001W	Rubber Tip	



Standard Steel Shower Stool

The basic steel shower stool is height adjustable. The plastic moulded seat is comfortable and hygienic.

Key features:

- Armrest assist with lowering and raising from seat
- Weight of stool 3.2kg
- AS/NZS 3973:2009
- ARTG 174421

Specifications:

Width between arms	45 cm (17¾")
Seat height	49 - 61 cm (19¼ - 24")
External depth	41 cm (16")
External width	51 cm (20")



*Now ordered in Box Qty

Part #	Description	Qty/Box
JAN-AC24-4	Standard Steel Shower Stool	4
JAN-AC20SP-004	Rubber Tip	1

Corner Shower Stool

The Corner Shower Stool features a handy tripod seat that fits into the corner of a shower with plastic moulded seat for comfort and hygiene. With a sturdy, height-adjustable frame and non-slip rubber feet, this aluminium stool is suitable for virtually any shower.

Key features:

- Sturdy, corrosion-resistant frame
- Height-adjustable legs with pin clips, and non-slip feet
- ARTG 174421

Specifications:

Seat width	52 cm (20½")
Seat depth	38 cm (15")
Height adjustment	44 - 54 cm (17¼ - 21¼")
No. of adjustment heights	5



Part #	Description
JAN-BE11	Corner Shower Stool

Aluminium Shower Stool with Padded Seat

The aluminium shower stool enables comfortable and safe showering. It is constructed from sturdy rust resistant aluminium with easy push button adjustment legs. The internal sliding legs provide smooth and effortless height adjustment whilst avoiding any chance of water filling the legs. Reinforced frame using 25mm tubing enables user to gain a firmer grip making transferring safer & easier.

Key features:

- Aluminium frame padded seat
- Height-adjustable legs with non-slip feet
- AS/NZS 3973:2009
- ARTG 291268

Specifications:

Seat width between arms	51.3 cm (20 ¹ / ₄ "
Armrest height from seat	13 cm (5")
Seat depth	25 cm (9 ³ / ₄ "
Seat height	52 - 67.7 cm (20 ¹ / ₂ - 26 ¹ / ₂ "
External depth	38.5 cm (15")
External width	53.3 cm (21")



*Now ordered in Box Qty

Part #	Description	Qty/Box
PAT-ACO6A-4	Shower Stool w Padded Seat	4



Extra Wide Heavy Duty Shower Stool with Padded Seat

An extra wide, comfortable shower stool with a padded seat. Corrosion resistant, zinc treated & epoxy powder coated steel frame.

Key features:

- Aluminium frame with padded seat
- Height-adjustable legs with non-slip feet
- Stacks vertically
- Weight of chair 3 kg (6¹/₂lbs)
- AS/NZS 3973:2009
- ARTG 291268

Specifications:

Seat width between arms	53.8 cm (21")
Armrest height from seat	14 cm (5 ¹ / ₂ "
Seat depth	27.5 cm (10 ³ / ₄ "
Seat height	53.9 - 68.9 cm (21 ¹ / ₄ - 27")
External depth	44.5 cm (17 ¹ / ₂ "
External width	65.5 cm (25 ³ / ₄ "
External height	67.9 - 82.9 cm (26 ¹ / ₂ - 32 ¹ / ₂ "



*Now ordered in Box Qty

Part #	Description	Qty/Box
PAT-ACO6-20HD-4	Extra Wide Heavy Duty Shower Stool w Padded Seat	4



Aluminium and Heavy Duty Shower Chair

Aluminium shower chair, which enables comfortable & safe showering, with a curved seat for comfort & drainage. Uncoated internal sliding legs provide smooth & effortless height adjustment, whilst avoiding any chance of water filling the legs. Reinforced frame using 25mm tubing enables users to gain a firmer grip making transferring safer and easier.

Key features:

- Lightweight and corrosion-resistant frame
- Stacks vertically
- Weight of chair 2.7 kg (6 lbs)
- AS/NZS 3973:2009
- ARTG 291268

Specifications:

Overall Height	77 - 92 cm (33¼ - 36¼")
Overall Width	59 cm (23")
Overall Depth	54.5 cm (25¼")
Seat Width	44 cm (17¼")
Seat Depth	40.5 cm (16")
Seat Height Range	41 - 56 cm (16¼ - 22")
Width between Armrests	51 cm (20")
Armrest (Seat to top of handle)	21 cm (8¼")
Backrest Height	35 cm (13¾")
Frame Diameter (armrest and leg)	25.4 mm (1")
Extension Leg Diameter	22.2 mm (¾"), HD 28.6 mm (1")
Rubber Tip	20.5 mm (¾"), HD 27 mm (1")

Standard

110 kg
17¼ st 242½ lbs
Maximum user weight

Heavy Duty

150 kg
23½ st 330½ lbs
Maximum user weight



*Now ordered in Box Qty

Part #	Description	Qty/Box
PAT-AC89A-3	Aluminium Shower Chair	3
PAT-AC89AHD-3	Aluminium Heavy Duty 150kg Shower Chair	3

Standard and Extra Wide Steel Shower Chair

The Steel Shower Chair is a standard sturdy chair, featuring a plastic contoured seat adds further comfort & assists drainage and is height adjustable. It has been cleverly designed with a gentle slope, so that water can flow away for easy drainage. The rubber feet are made from non slip rubber. Manufactured from corrosion resistant, zinc treated and epoxy powder coated steel. Assembled and stacks vertically.

Key features:

- Contoured plastic moulded seat
- Stacks vertically
- AS/NZS 3973:2009
- ARTG 291268

Specifications:

Overall Height	77 cm - 92 cm (33¼ - 36¼")
Seat Width	44 cm (17¼")
Seat Depth	40.5 cm (16")
Seat Height Range	42 cm - 57 cm (16½ - 22½")
Armrest (Seat to top of handle)	21 cm (8¼")
Backrest Height	35 cm (13¾")
Frame Diameter (armrest and leg)	22.2 mm (¾")
Extension Leg Diameter	25.8 mm (1")

Standard

Extra Wide

Overall Width	59 cm (23")	64 cm (25¼")
Overall Depth	54.5 cm (25¼")	59.6 cm (23½")
Width between Armrests	51 cm (20")	58 cm (22¾")
Weight	3.7kg (8.2 lbs)	4.7kg (10.4 lbs)

125 kg
19 st 264½ lbs
Maximum user weight



*Now ordered in Box Qty

Part #	Description	Qty/Box
PAT-70495-3	Standard Steel Shower Chair	3
PAT-AC89-20-3	Extra Wide Steel Shower Chair	3





Aluminium Bedside Commode Chair

The Aluminium Bedside Commode is comfortable and hygienic yet modest. Features rust resistant aluminium frame, which is sturdy and height adjustable. This commode chair features non slip rubber tips and an easily removable backrest. The lift out bowl with lid, padded vinyl contoured backrest, toilet seat and cover offer extra comfort and pressure relief whilst maintaining ease of cleaning & hygiene.

Key features:

- Lightweight, corrosion-resistant aluminium frame
- Bowl and lid included
- 12 month warranty
- AS/NZS 3973:2009

Specifications:

Seat width between arms	45 cm (17¾")
Overall depth and width	55 cm (21")
Seat height	44 - 59 cm (17¼ - 23¼")
Armrest height from seat	18 cm (7")



Bedside Commode Chair

The Bedside Commode has a clip on toilet seat, which clips on in seconds. There are no metal fittings on toilet seat. It is also height adjustable and easy to clean. It has a light blue deluxe padded overseat, as well as a removable round bowl with lid & handle. The lid eliminates spillages and reduces risk of cross infection. This bedside commode looks very smart in a bedroom.

Key features:

- Powdercoated dark blue frame
- Padded seat and backrest
- Bowl and lid included
- Stackable up to x10
- 12 month warranty

Specifications:

Seat height	46 - 61 cm (18 - 24")
External height	82 - 97 cm (32¼ - 38¼")
External length	55 cm (21½")
External width	61 cm (24")



*Now ordered in Box Qty

Part #	Description	Qty/Box
NOV-COM71137-4	Aluminium Bedside Commode	4
MFI-AC90B	Bowl with Lid - Grey	
NOV-SP-V72002	Padded Toilet Seat	
NOV-SP-N8590-B11	Replacement Arm Rests - Pair	
NOV-CHA71237	Replacement Backrest	
NOV-BAT70373	Rubber Tips	
NOV-MOBWAL8590	Screw Pack	

*Now ordered in Box Qty

Part #	Description	Qty/Box
SHH-DAYCOM521A-3	Bedside Commode	3
SHH-521APL002	Ringseat - Black	
SHH-521ASA003	Padded Seat	
MFI-AC90B	Bowl with Lid - Grey	
SHH-521A-000003	E Clip - Grey	
SHH-521A-000004	Rubber Sleeve - Grey	
SHH-521A-000005	Rubber Tip - Grey	



Deluxe Shower Commode Chair

The Deluxe Shower Commode Chair features a lightweight, aluminium frame, making it rust-free and highly durable. Attendant- or self-propelled, this chair offers independence for bathing and toileting. With versatile swing-away, detachable footrests and flip-up arms, it is ideal for use in a shower cubicle or wet room, as a commode or over a standard toilet bowl. The detachable back enables the chair to be partially dismantled for more compact storage.

Key features:

- Lightweight, aluminium frame with cut-away seat
- Seat insert
- Brakes on all four wheels
- Includes square commode pan, footrests and seat insert
- ARTG 199705

Specifications:

	Attendant-propelled	Self-propelled
Depth without footrests	55 cm (21¾")	84 cm (33")
Depth with footrests	81 cm (31¾")	111 cm (43¾")
Overall width	55.5 cm (21¾")	69 cm (27")
Rear self-propelling wheels		60 cm (24")
Weight	12.4 kg (27.3 lbs)	20.7 kg (45.6 lbs)
Seat height	49.5 - 57 cm (19½ - 22½ ")	
Seat width	55.5 cm (21½")	
Seat depth	42.5 cm (16¾")	
Overall height	95.8 - 100.8 cm (37¾ - 39¾")	
Aperture size	19 x 36 cm (7½ x 14")	
Front wheel diameter	13 cm (5")	
Under seat clearance	43 - 50.5 cm (17 - 19")	



Adjustable Height Attendant-propelled

NEW-091563451



Standard Self-propelled

NEW-081135367

Part #	Description	Seat width
NEW-091563451	Attendant-propelled, adjustable height	54 cm (21¼")
NEW-091183698	Self-propelled, fixed height	54 cm (21¼")

Deluxe Tilt n Space Shower Commode Chair

The Deluxe Tilt n Space Shower Commode Chair includes all the great features of the Deluxe Shower Commode Chair but has the additional advantages of being height-adjustable and tilting to provide a secure and comfortable position for patients with more complex support requirements. It is a sturdy commode that can tilt the user to provide a secure and comfortable position for washing and toileting. It can be reclined from -5° to 30° with padded neck support that can be positioned at any angle for maximum comfort.

Key features:

- Lightweight, aluminium frame with cut-away seat and seat insert
- Swing-away, detachable footrests and flip-up arms
- Brakes on all four wheels
- Includes square commode pan, footrests, seat insert and adjustable neck support

Specifications:

Seat width	45.7 cm (18")
Seat depth	42.5 cm (16¾")
Seat height	53 - 60.5 cm (20¾ - 23¾")
Overall width	56 cm (22")
Overall height	117 cm (46")
Aperture size	19 x 36 cm (7½ x 14")
Wheel diameter	13 cm (5")
Under seat clearance	43 - 50.5 cm (17 - 20")
Weight	10 kg (22 lbs)



Part #	Description
NEW-091563501	Deluxe Tilt n Space Shower Commode Chair

Deluxe Heavy Duty Shower Commode Chair


The Deluxe Shower Commode Chair features a lightweight, aluminium frame, making it rust-free and highly durable. This chair offers independence for bathing and toileting. With versatile swing-away, detachable footrests and flip-up arms, it is ideal for use in a shower cubicle or wet room, as a commode or over a standard toilet bowl. The detachable back enables the chair to be partially dismantled for more compact storage.

Key features:

- Lightweight, aluminium frame with cut-away seat and seat insert
- Heavy-duty version features increased seat width
- Brakes on all four wheels
- Includes square commode pan, footrests and seat insert

Specifications:

Seat height	49.5 - 57 cm (19½ - 22½")
Seat depth	42.5 cm (16¾")
Overall height	99 cm (40")
Depth without footrests	55 cm (21¾")
Depth with footrests	81 cm (31¾")
Aperture size	19 x 36 cm (7½ x 14")
Front wheel diameter	13 cm (5")
Rear self-propelling wheels	60 cm (24")
Under seat clearance	43 - 50.5 cm (17 - 20")



Part #	Description	Seat width
NEW-091563477	Attendant-propelled	56 cm (22")



Tilt n Space Attendant-propelled

NEW-091563501



Heavy-duty Attendant-propelled

NEW-091563477



Long-handled Sponge

The Long-handled Sponge is mounted on a plastic-coated wire to help reach the legs, feet and back. It is ideal for patients with limited reach.

Key features:

- Plastic-coated, flexible wire
- Oval foam sponge
- Grey plastic built-up handle



Part #	Description	Weight
PAT-AA1831CY	Length 61 cm (24")	100 g (3.5 oz)
PAT-AA1831B	Length 38 cm (15")	75 g (2.75 oz)

Long-handled Sheepskin Pad

A luxury double thickness sheepskin pad with a plastic built-up handgrip. Can be used for washing or powdering. Pad size 95 x 125mm.

Key features:

- Soft, sheepskin pad
- Plastic-coated, flexible wire
- Grey plastic built-up handle



Part #	Description	Weight
PAT-AA1841C	Length 61 cm (24")	135 g (4.8 oz)

Long-handled Toe Washer

The Long-handled Toe Washer has a plastic built-up handle and a flat wire frame covered in a white towelling pad for washing between the toes.

Key features:

- Plastic-coated, flexible wire
- Grey plastic built-up handle
- Two washing pads included

Part #	Description	Weight
PAT-AA1847	Length 71 cm (28")	120 g (4.2 oz)
PAT-AA184601	Replacement Pads - Pair	



Inflatable Shampoo Basin

The Inflatable Shampoo Basin is designed for comfortable in-bed shampooing with secure, form-fitting support that cushions the head and shoulders. It is lightweight, portable and compact for convenient storage.


Key features:

- Built-in inflatable headrest
- Made from heavy-duty vinyl
- Flexible drainage hose includes stopper for control when emptying

Specifications:

Width	61 cm (24")
Depth	20 cm (8")
Length	51 cm (20")
Weight	0.68 kg (1½ lbs)

Part #	Description
PAT-AA1874	Inflatable Shampoo Basin



Flannel Bath Strap and Wash Mitt

The dual-sided Flannel Bath Strap and Wash Mitt helps the bather reach awkward places. Use the cotton terry surface to wash and the coarse foam reverse side to scour and tone skin. The D-shaped handles create an easy, secure grip, making washing easier.

Key features:

- Cotton flannel bath strap with D-shaped plastic handles
- Set includes matching hand mitt
- Hand-washable
- Strap dimensions: 80 x 10 cm (31 ½ x 4"), total length (including handles) 94 cm (37")

Part #	Description	Weight
PAT-AA1856	Flannel Bath Strap and Wash Mitt	100 g (3.5 oz)



Dual Function Lotion Applicator

The long-handled Dual Function Lotion Applicator allows users to apply cream and massage their skin simultaneously. Rotating balls ensure even distribution and a pleasant massage. An interchangeable foam head helps when applying thicker creams directly onto the skin. A long, moulded handle makes reaching the back and legs easier. Great for sunscreen application.

Key features:

- Dual-sided applicator head
- Moulded plastic handle
- No replaceable heads

Part #	Description	Weight
PAT-AA1818	Handle length 35 cm (13¾")	150 g (5.3 oz)



Lotion Applicator

Specially designed to help apply moisturising lotions, suntan oils, sport creams and medications on hard to reach places. The built up, ribbed handle has a thumb rest and is angled in two places for maximum manoeuvrability. The sponge can be removed when worn and replaced with a new one.

Key features:

- Ribbed handle with thumb rest
- Removable sponge with replacement available

Part #	Description	Weight
PAT-AA1862	Handle Length 30.5 cm (12")	50 g (1.8 oz)
PAT-AA186201	Replacement Sponge - Pair	



EasyBar Suction Grab Rail

The EasyBar Suction Grab Rail offers sturdy and safe support to the user without having to be permanently attached to the wall. The rails assist individuals with sit-to-stand transfers, helping increase safety, balance and stability. The non-slip handle offers greater independence by providing stability when getting in and out of the bath or shower or for general use around the home. Safety lock indicators display green when the rail is securely attached and red when unsecured. The EasyBar was designed for use only on smooth and nonporous surfaces such as tiles, glass, porcelain or fibreglass.

Key features:

- Equipped with two pressure switches on each end that create a strong vacuum to latch onto surfaces
- Suitable for temporary or permanent use
- Available in two lengths

Part #	Description
PAT-091325786	35 cm (13¾")
GAL-ADLBAT75680	50 cm (19½")



Tub Mount Bath Rail

The Bath Grab Rail provides a firm handle at the bath side, helping to give confidence and stability when entering and exiting the bath. The coated steel rail is fixed to the side of the bath by a clamp and cushioned on the inside by rubber pads.

Key features:

- Coated steel frame fits most baths
- Screw-tight clamp with rubber pads
- Not for use with lipped-style baths

Specifications:

Length	45.7 cm (18")
Height above bath side	34.9 cm (13¾")
Clamp width adjustment	5.5 - 11.5 cm (2¼ - 4½")
Weight	2.4 kg (5¼ lbs)

100 kg

15¾ st220½ lbs

Maximum user weight

*Now ordered in Box Qty

Part #	Description	Qty/Box
JAN-JMC-5420-4	Tub Mount Bath Rail	4



Savanah Raised Toilet Seat

The Savanah Raised Toilet Seat features a gently contoured surface for extra comfort and slopes down slightly from the back to the front for easier lowering and raising. A large rear cut-out gives clearance for the coccyx and combined with the gentle downward slope of the seat, fosters rear personal hygiene and an expansive front cut-out. The deep, built-in splash guard minimises soiling. Fully sealed, the seat is resistant to stains and odours and is easy to clean.

Two L-shaped moulded plastic brackets have bonded non-slip rubber buffers to keep the toilet firmly in place. The brackets are easily adjusted to fit most toilet shapes and sizes using two large, fluted turning knobs, making fitting and removing the seat quick and easy. No front bracket or front-turning knob to avoid snagging skin or clothes.



Key features:

- Lightweight, plastic construction for strength and durability
- Fully sealed, resistant to stains and odours
- Available in three seat heights
- Available with or without lid

Specifications:

Overall length	40.7 cm (16")
Overall width with brackets	37 cm (14¾")
Aperture	26.6 x 20.7 cm (10½ x 8¼")



5cm



10cm



15cm

Part #	Description	Height	Weight
PAT-AA2112	Without Lid	5 cm (2")	0.8 kg (1¾ lbs)
PAT-AA2114Y	Without Lid	10 cm (4")	1 kg (2¼ lbs)
PAT-AA2116	Without Lid	15 cm (6")	1.15 kg (2½ lbs)
PAT-AA2112L	With Lid	5 cm (2")	1.2 kg (2¾ lbs)
PAT-AA2114L	With Lid	10 cm (4")	1.4 kg (3 lbs)
PAT-AA2116L	With Lid	15 cm (6")	2 kg (4½ lbs)



Raised Toilet Seat with Armrests

Best seller, feature packed, improved design. The Raised Toilet Seat with Armrests fosters increased independence, notably by compensating for impaired movements. The seat with detachable armrests comes with a lid that can be removed for easy cleaning. The versatile design includes a choice of three fixed seat heights to suit the user’s needs. Arms swing back to aid transfer on and off the seat.

Key features:

- Moulded plastic base includes lift-up ABS armrests and a steel frame
- Cross bars adjust in width for each seat position
- Existing toilet seat must be removed prior to installation

Specifications:

Seat width	36 cm (14")
Seat depth	42.5 cm (16⅔")
Overall width	50 cm (19½")
Overall depth	47.5 cm 18⅔")
Overall height	38 cm (15")
Seat height adjustments	6 cm (2"), 10 cm (4") and 15 cm (6")
Aperture size	21 x 27.5 cm (8¼ x 10⅔")
Weight	4.3 kg (9½ lbs)

150kg

23½ st 330 lbs

Maximum user weight



Part #	Description
PAT-091558873	Raised Toilet Seat with Armrests

Height-adjustable Toilet Surround

The Height-adjustable Toilet Surround, in steel provides extra support and safety when using the toilet. The height-adjustable frame features moulded plastic armrests for additional comfort.

Key features:

- Available with steel frame
- Height-adjustable legs with non-slip feet
- Moulded plastic armrests

Specifications:

Height	69.5 - 85 cm (27½ - 22½")
Width	59 cm (23¼")
Depth	46 cm (18")
Weight	3 kg (6¾ lbs)

160 kg

25 st 352¾ lbs

Maximum user weight

*Now ordered in Box Qty

Part #	Description	Qty/ box #
HOM-503A-5	Height-adjustable Toilet Surround — Steel	5
HOM-SP-DAYHOMCCLIP	C-Clip for Toilet Surround	



Devon Folding Support Rail

The Devon Folding Support Rail is a sturdy folding rail that projects at right angles to the wall. The rail is held in the upright position by a spring-loaded cam and is easily folded down.

Key features:

- White, powder-coated steel rail
- Attaches to wall with steel wall bracket
- Available in standard and short length
- ARTG 210113

Specifications:

Rail diameter	2.5 cm (1")
Backplate dimension	30 x 12.5 cm (11¾ x 5")

67 kg
10½ st 147¼ lbs
Maximum load weight

Part #	Rail length	Weight
PAT-AA2002	Standard — 76 cm (30")	3.0 kg (6½ lbs)
PAT-AA2000	Short — 55 cm (21½")	2.6 kg (5¾ lbs)



Devon Deluxe Folding Support Rail

The Devon Deluxe Folding Support Rail offers the same features as the Devon Folding Support Rail but with a larger diameter rail for easier gripping, a spring mechanism to hold it upright, as well as a locking grip. Available as a standard folding rail with a supporting leg or as a floor-mounted rail, this product is available with multiple options to meet the user's specific requirements.

Key features:

- White, powder-coated steel rail with larger diameter rail
- Spring mechanism and lock hold rail in the upright position
- ARTG 210113

Specifications:

Rail diameter	3.2 cm (1¼")
Backplate dimension	30 x 12.7 cm (11¾ x 5")

84 kg
13¼ st 185¼ lbs
Maximum load weight

Part #	Rail length	Weight
PAT-AA2012	Standard — 77 cm (30⅓")	4.1 kg (9 lbs)
PAT-AA2028	Toilet Roll Holder	110g (3.9 oz)



Devon Deluxe Floor-mounted Folding Rail

The Devon Deluxe Floor-mounted Folding Rail has a large diameter rail for easier gripping and a spring mechanism to hold it in the up position. It is mounted on a floor-fixing post for positioning versatility with a clamp to adjust by height.

- ARTG 210113

Specifications:

Floor plate	25 x 15 cm (10 x 6")
Floor post height	111 cm (44")
Post height adjustment	58.4 - 101.2 cm (23 - 40")

127 kg
20 st 280 lbs
Maximum load weight

Part #	Rail length	Weight
PAT-AA2026	Standard — 76.5 cm (30")	8.4 kg (18½ lbs)



Folding Over Toilet Aid

The Folding Over Toilet Aid assists transferring on & off the toilet. The folding frame ensures compactness in small spaces. Constructed of rust resistant aluminium, the frame is height adjustable & weighs just 3.8kg (one of the lightest available). The non slip rubber feet offer stability & steadiness & the armrests offer additional support for standing. The plastic seat is hygienic & easy to clean. The Folding Over Toilet Aid is an ideal multi-purpose piece of equipment. Perfect for travelling, for use as a shower chair or used with bowl & lid as a bedside commode.

- ARTG 199500

Specifications:

Seat width between arms	48 cm (18¾")
Seat depth	41 cm (16¼")
External width	53.5 cm (21")
Overall seat depth	50.2 cm (19¾")
Armrest height from seat	22 cm (8½")
Height clearance	42 - 58 cm (16½ - 22¾")
Width clearance	42 cm (16½")
Weight	3.65 kg (8 lbs)



Part #	Description
NOV-8560-A	Folding Over Toilet Aid (without bowl and lid)
NOV-BAT68662	Bowl and Lid
NOV-8560-B01	Splashguard
NOV-BAT70430	Toilet Seat

Part #	Description
NOV-BAT71089	Rubber Tips
NOV-BAT72012	Rubber Seat Stoppers
NOV-BAT70378	Screw Set (4 screws and washers)

Aluminium Over Toilet Frame

For those who have difficulty bending or lowering themselves, particularly after a hip replacement, raising the height of the toilet seat makes toileting much safer & easier. Over Toilet Frame with a seat & lid, have arms for added stability, support & grip, while the legs are height adjustable. Constructed from sturdy rust resistant aluminium. The plastic seat & lid clip off to facilitate hygienic & easy cleaning. The anti-slip rubber tips offer stability & steadiness even on wet surfaces. Easy push button height adjustment legs. The uncoated internal sliding legs provide smooth & effortless adjustment whilst avoiding any chance of water leakage into the legs. The reinforced frame uses 25mm tubing, enabling users to gain a firmer grip whilst making transfers.

ARTG 199500

Specifications:

Seat width between arms	45 cm (17¾")
Seat height range	44 - 59 cm (17¼ - 23¼")
External width	55 cm (21½")
External depth	55 cm (21½")
Armrest height from seat	18 cm (7")
Height clearance around toilet	48 cm (18¾")
Height clearance over toilet	40 - 55 cm (15¾ - 21½")
Weight	2.7 kg (6 lbs) : with lid 3.4 kg (7½ lbs)



*Now ordered in Box Qty

Part #	Description	Qty/box #
NOV-AC02A-6	Over Toilet Frame	6
NOV-BAT68662	Bowl and Lid	
NOV-8560-B01	Splashguard	

Part #	Description
NOV-N8500CA17	Toilet Seat
NOV-N8500CA18	Toilet Lid
NOV-BAT70373	Rubber Tips
NOV-BAT72012	Rubber Seat Stoppers

Stirling Width Adjustable Toilet Frame

The Stirling Toilet Frame features a sturdy, width adjustable powder-coated metal frame with a clip-on, one-piece moulded toilet seat that is slightly angled forward to make lowering and raising easier. The contoured seat has a partial front cut-away to facilitate personal hygiene. The width-adjustable option features a pin-clip system, enabling the frame to fit a variety of toilet styles and be completely assembled or disassembled without tools. The deluxe version includes padded armrests on the side rails for additional comfort and easier grip.

Key features:

- Frame legs feature non-slip rubber feet
- Pin-clip system for height and width adjustments leaves no protrusions to snag skin or clothing
- Flat-packed, assembly required
- Available in three versions: standard, width-adjustable or width-adjustable with padded arms

165 kg
26 st 363½ lbs
Maximum load weight



Part #	Description
PAT-AA2214	Standard Height-/Width-adjustable Frame
PAT-AA2218	Floor Fixing Feet - Set of 4
PAT-AA2251	Delux Height-/Width-adjustable Frame with Padded Arms

Over Toilet Aid

The Days Over Toilet Aid assists transfers on and off the toilet. This toilet surround frame is height adjustable, with a push button adjustment. Powder coated steel frame with non slip rubber feet. Comes complete with a seat and splashguard.

Specifications:

Seat width	37 cm (14½")
Width between arms	43 - 57 cm (16¾ - 22½")
Height adjustable	44 - 59 cm (17¼ - 23¼")

125 kg
19½ st 275½ lbs
Maximum user weight

*Now ordered in Box Qty

Part #	Description	Qty/box #
JAN-AC20-4	Days Over Toilet Aid	4
JAN-AC20SP-008	Toilet Seat	
JAN-AC20SP-009	Bowl / Splashguard	
JAN-AC20SP-004	Rubber Tip	



Bariatric All in One

Redesigned with the carer in mind, the Bariatric All In One is an over toilet aid, shower chair and commode. Comes with square bowl with easy front access to remove bowl. The square bowl does not have a lid. Optional splash guard available. The seat lifts up and comes with a seat cover. The removable backrest adds stability when used as a bedside commode and can be removed when using over a toilet. The moulded armrests assist gripping for safe and easy transferring.

Specifications:

Width between arms	60 cm (23½")
Height adjustable	44 - 59 cm (17¼ - 22½")

295 kg
46½ st 650 lbs
Maximum user weight

*Now ordered in Box Qty

Part #	Description	Qty/box #
MAX-AUSBAT70282-2	Bariatric All in One	2
MAX-70282-01	Toilet Seat with Hole	
MAX-70282-02	Toilet Seat Cover without Hole	
MAX-70282-03	Square Bowl - No Lid	
MAX-70282-05	Rubber Tips	



Sherwood Bariatric Commode

The Sherwood Bariatric Commode offers the user comfort, stability and independence for use anywhere within the home. The durable metal frame is height-adjustable and designed for bariatric users, featuring a thick padded backrest, toilet seat and seat cover, and removable armrests that allow for easy transfers from the bed, wheelchair or armchair.

Key features:

- Coated metal frame with additional crossbar support
- Padded toilet seat and backrest
- Height-adjustable legs include non-slip feet
- Includes square commode pan and separate seat cover

Specifications:

Seat width	66 cm (26")
Seat depth	50 cm (20")
Seat height	48 - 63.5 cm (19 - 25")
Height of armrest	22 cm (8¾")
Size of aperture	28 x 23 cm (11 x 9")
Footprint	67 x 59 cm (26½ x 23¼")

254 kg
40 st 560 lbs
Maximum user weight

Part #	Description
PAT-091422963	Sherwood Bariatric Commode



Folding Bottom Wiper

The Folding Bottom Wiper features a moulded plastic head with recessed serrations to grip paper. It folds into three sections for easier storage. Designed to be used from the front.

Specifications:

Length	26.5 cm (10½")
Fold-down length	12 cm (4¾")
Weight	0.8 kg (1¾ lbs)



Part #	Description
PAT-AA2652	Folding Bottom Wiper

Bottom Wiper

Designed by the Bath Institute for Medical Engineering, the Bottom Wiper has a moulded plastic head with recessed serrations to grip paper. It is mounted on a long, lightweight, curved handle that assists in reaching areas where wiping is difficult.

Specifications:

Length	26.5 cm (10½")
Weight	0.7 kg (1½ lbs)



Part #	Description
PAT-AA2650	Bottom Wiper

Savanah Splash Guard

Simple to fit, the Savanah Splash Guard clips onto a toilet bowl, a raised toilet seat or a toilet frame to prevent embarrassing and unhygienic accidents. Simple to remove for cleaning and getting on and off the toilet.



Part #	Description
PAT-AA2602	Savanah Splash Guard

General Purpose Bedpan

The General Purpose Bedpan is made from polypropylene and can be cleaned with steam, hot water or chemicals and has a 1-litre capacity.

Specifications:

Width	30 cm (11¾")
Length	38 cm (15")
Depth	10 cm (4")
Weight	1.1 kg (2½ lbs)



Part #	Description
PAT-091081876	General Purpose Bedpan

Economy Male or Female Urinal

The Economy Male Urinal is moulded from transparent plastic and features a snap-on lid and handle. The long neck ensures it is easy to use. It has a 1-litre capacity.



PAT-AA2668

The Economy Female Urinal has a handle for ease of use. It is anatomically designed to use when sitting, standing or lying down. It has an 800-mL capacity with graduations for output.



PAT-AA2666

Part #	Description	Weight
PAT-AA2668	Male Urinal	0.1 kg (3.5 oz)
PAT-AA2666	Female Urinal	0.1 kg (3.5 oz)

Inflatable Rubber Ring Cushion

Designed for those who are temporarily or permanently confined to bed or chair and are thus prone to bedsores, discomfort or pain. Also useful for those with sedentary occupations where conditions such as haemorrhoids and varicose veins can arise from inhibited circulation.



Part #	Description
PAT-AA3938	Inflatable Rubber Ring Cushion





Aids to Daily Living



Dressing Stick

The Dressing Stick helps individuals put on and take off garments they cannot easily reach or that require movements that are difficult, such as socks, jumpers and jackets. The Standard dressing stick is equipped with a rubber tip on one end and a double-wire hook on the other. The Deluxe dressing stick has a large dressing hook at one end and a smaller looped hook at the other end.

Specifications:

Length	52 cm (20½")
Weight	36 g (1.3 oz)

Part #	Description
PAT-AA4680	Dressing Stick : 52 cm (20½")
PAT-A6983	Deluxe Dressing Stick : 65 cm (25½")



Kings Button Hook

The Kings Button Hook is a shaped stainless-steel wire constructed to assist with one-handed dressing. Simply pass the hook through a buttonhole and over the button, then pull the button back through the buttonhole.

Specifications:

Wire length	9.7 cm (3¾")
Handle	12 cm (4¾")
Weight	75 g (2.6 oz)



Part #	Description
PAT-AA4677	Kings Button Hook

Button Hook and Zip Puller

The Button Hook and Zip Puller is a 2-in-1 tool that assists individuals in buttoning clothing and pulling zippers. A button hook on one end and a zip puller on the other help with various dressing tasks.

Specifications:

Length	24.5 cm (9¾")
Handle width	1.5 cm (¾")
Wood length	17 cm (6¾")
Weight	8 g (0.3 oz)



Part #	Description
PAT-AA4673	Button Hook and Zip Puller




Plastic Shoehorn

The Plastic Shoehorn features a long, curved handle designed to extend reach when putting on shoes. Includes an attached loop for storage.

Specifications:

Length 43 cm (17")
Weight 37 g (1.3 oz)



Rigid Leg Lifter

The Rigid Leg Lifter is a webbing strap with a stiffened stem that enables the user to move an immobile leg when in a bed or wheelchair or using a footstool. The upper loop fits around the hand and wrist. The lower loop slides over the foot to offer control in lifting the leg.

Specifications:

Length 86.4 cm (34")
Weight 78 g (2.7 oz)



Part #	Description
PAT-091095157	Plastic Shoehorn

Part #	Description
PAT-AA8620	Rigid Leg Lifter

Chrome-plated Metal Shoehorn

The Chrome-plated Metal Shoehorn features a chrome metal finish on a long-stemmed shoehorn. The blade is mounted on a spring to allow it to bend in relation to the handle when putting on shoes. The plastic-moulded handle has a comfortable grip, making it easy for the user to hold.

Specifications:

Length 60 cm (23¾")
Weight 115 g (4 oz)



Dual-Handle Leg Lifter

The Dual Handle Leg Lifter is designed to help users reposition their legs when they have limited or no use of their legs. Two handles enable users to utilise both arms when moving the leg, and the sturdy centre rod stabilises the device. Ideal for using while in a chair or bed.

Specifications:

Length 86.5 cm (34")
Weight 110 g (3.9 oz)



Part #	Description
PAT-091095199	Chrome-plated Metal Shoehorn

Part #	Description
PAT-AA8620B	Dual-Handle Leg Lifter

Coated Metal Shoehorn

The Coated Metal Shoehorn features a sturdy, long design with epoxy-coated steel and a soft, hook-shaped PVC handle.

Specifications:

Length 60 cm (23¾")
Weight 233 g (8.2 oz)



Coiler Shoe Laces

These elastic based laces do not require tying. Lace the shoes up and release the laces at the required pressure and the lace will coil to keep in place. Ideal for those who find it painful to bend and tie traditional laces.



Part #	Description
PAT-AA4606Y	Coated Metal Shoehorn

Part #	Description	Colour
PAT-AA4618A	Coiler Shoe Laces	White Pair
PAT-AA4618C	Coiler Shoe Laces	Brown Pair
PAT-AA4618G	Coiler Shoe Laces	Black Pair



Terry Cloth Sock Aid

The Terry Cloth Sock Aid makes sock application easy. The flexible, plastic sock aid design enables the aid to bend around the heel, making it easier for the user to pull socks on over their feet. Nylon lines the inside to reduce friction, and terry cloth lines the outside to hold the sock firmly in place. Two pulling tabs with loop handles are ideal for one-handed use.

Specifications:

Length	23 cm (9")
Weight	71 g (2.5 oz)




Part #	Description
PAT-A7273	Terry Cloth Sock Aid

Brevetti Sock and Stocking Aid

The Brevetti Sock and Stocking Aid is designed to help those with restricted mobility to apply socks and stockings easily. The Brevetti consists of a coated steel tube with a plastic collar to hold the garment open on one end and a plastic shoehorn with a notch for easy sock removal on the other.

Specifications:

Length	73.7 cm (29")
Weight	386 g (13.6 oz)




Part #	Description
PAT-AA4666	Brevetti Sock and Stocking Aid

Dorking Stocking Aid — Standard or Side Slot

The Dorking Stocking Aid is a simple, inexpensive device that assists the user when putting on stockings or socks, especially when bending to reach the foot is difficult. It features a flexible plastic gutter with long cotton tapes attached to the top. Feed the stocking or sock through the plastic and hold it in place with the side notches. Place the foot in the open end of the stocking and pull it over the foot and up the leg using the tapes. Available with or without side slots.

Specifications:

Cord length	92 cm (36")
Gutter length	22.9 cm (9")
Weight	55 g (2 oz)



Part #	Description
PAT-AA4650Y	Standard
PAT-AA4652Y	Side Slot

Dorking Stocking Donner

The Dorking Stocking Donner is a compact, durable tool that extends the user's reach to help them dress independently. The extendable, hold-down handles rotate 90°, creating a comfortable position for the user while applying stockings. The central column, which holds the stocking, has a small ridge that prevents movement and widens at the top to accommodate the larger calf. The epoxy coating creates a smooth surface, reducing the risk of snagging.

Specifications:

Height	31 cm (12¼")
Width	20.5 cm (8")
Column height	20.5 cm (8")
Column width	12.5 cm (5")
Column depth	8.5 cm (3½")
Height with handles	56 cm (22")
Weight	1.1 kg (2½ lbs)



Part #	Description
PAT-AA4653	Dorking Stocking Donner

Compression Stocking Aid

A stocking aid which will apply support hosiery quickly and evenly. Grooves catch the compression stocking and prevent it from popping off. Works like a traditional sock aid - just pull the stocking over the shell, insert the foot and pull the cord to bring the hose up onto the leg. Suitable for Class 1 and 2 stockings. Recommended for below knee stockings.



Specifications:

Length	14.5 cm (5¾")
Width	15 cm (6")

Part #	Description
PAT-AA75441	Sammons Preston Compression Stocking Aid

Step Stool

The Step Stool is a simple, sturdy stool consisting of a chrome-plated steel frame with non-slip rubber mats. Rubber ferrules provide stability on floor surfaces. Available with or without a rail.



Specifications:

Top area	28 x 36 cm (11 x 14¼")
Base area	30.5 x 41 cm (12 x 16¼")
Base height	22.5 cm (9")
Maximum load	120 kg (18¾ st; 264½ lbs)

Part #	Description
JAN-JMC-601	Step Stool
JAN-JMC-5710	Step Stool with Handrail

Easi-Grip Scissors

Easi-grip scissors are extremely lightweight with plastic handles and stainless steel blades. The continuous loop handle automatically reopens the scissors when pressure is released. As they require only a gentle squeeze to operate them, they can be used between fingers and thumb or fingers and the palm of the hand. Available with round or pointed blades. They can be used either right or left-handed.




Part #	Blade Length	Style	Weight
PAT-AA6710	45 cm (17¾")	Round End	30g
PAT-AA6714	45 cm (17¾")	Pointed	30g
PAT-AA6716	75 cm (29½")	Pointed	46g

Dorking Tights Aid

This works on the same principle as the stocking aid, but has a double gutter to enable it to be used for tights.

Specifications:

Cord Length	92 cm (36¼")
Length	22.9 cm (9")
Weight	106 g (3.7 oz)



Part #	Description
PAT-AA4654Y	Dorking Tights Aid

Long Handled Combs and Brush

These utensils assist those who cannot raise their arm or reach behind their head to comb their hair. They consist of a lightweight moulded plastic handle with a push-on extension piece, a hinged end piece which holds the comb or brush and has two different angled positions. When not in use, the comb or brush folds back into the handle and the extension can be taken off to give a pocket sized item.

Specifications:

Length folded	23.5 cm (9¼")
Length with handle	32 cm (12½")
Length with extension	46 cm (18")



Part #	Description	Weight
PAT-AA4712	Style Comb	114 g
PAT-AA4714	Brush	128 g

Easi-Grip Chiropodist Scissors

The long loop handles allow the whole hand to be used for controlling the scissors. This is ideal for those with a weak grip. The extended shank helps users who have difficulty bending. The angled blades give a good cutting position, with the top blade being slightly serrated to stop it slipping. Stainless steel blades. Overall length 200mm. Blade length 30mm. Weight 55g.



Part #	Description
PAT-AA6766Y	Easi-Grip Chiropodist Scissors



Key Turners

The Key Turners have generous curved built up handles for holding keys. They provide an easy grip and good leverage for turning any Yale or mortice type door key. The two key turner feature a locking lever that allows keys to be folded into the handle when not in use.



Part #	Description	Colour
PAT-AA6242Y	Key Turner Double	Blue

Medi-Dispenser

The Medi-Dispenser is an organizer that holds a week's medication with as many as four daily doses. The plastic wallet has seven plastic containers printed with the day of the week and embossed with the number of the day of the week.



Part #	Description
PAT-AA9900Y	Medi-Dispenser

Contour Turner

The Contour Turner offers leverage for turning small household items, such as gas or radiator taps. The head contains spring-loaded stainless-steel rods, which when pressed around an object, will retract and conform around it, creating a grip around the object when the handle is turned. The handle is made of easy-to-grip moulded plastic.



Specifications:

Head size — width	9.5 cm (3¾")
Head size — depth	5 cm (2")
Weight	130 g (4.6 oz)

Part #	Description
PAT-AA6230	Contour Turner

Pen and Pencil Holder

The Pen and Pencil Holder is constructed from soft PVC moulding, which slides onto pens and pencils to make gripping and controlling writing easier. Supplied in a pack of three.




Part #	Description
PAT-AA7200SY	Pen and Pencil Holder - Set of 3

Plug Pull

The Plug Pull is a self-adhesive, plastic moulded T-piece that fixes to the back of a standard, three-pin plug, providing leverage to remove the plug from the socket. It is also helpful for drawers and doors where extra leverage is required. Supplied in a pack of three.



Part #	Description
PAT-AA6310	Plug Pull - Set of 3



Pickup Reacher

The Pickup Reach is a lightweight, robust reacher ideal for retrieving or reaching various items in the home, garden or stores. Available in four lengths and with optional extras, it has twin-gripping jaws with sculptured liners and a trigger designed for comfort. A magnetic head enables easy retrieval of pins, paper clips and other small, metal items. The handle is large and roomy to allow for the use of all fingers, and the short trigger distance is ideal for those with stiff, arthritic fingers or limited mobility. It has an easily rotatable head and a removable stick clip to connect the reacher to wheelchairs, walking frames and other mobility devices.

Specifications:

	Length	Weight
Standard	60 cm (24")	215 g (7.6 oz)
Long	75 cm (30")	235 g (8.3 oz)
Extra long	90 cm (35½")	250 g (8.8 oz)



Trigger Design



The trigger is designed for comfort and is large and roomy to allow use of all fingers. The short trigger distance is ideal for stiff, arthritic fingers. The trigger projects through the top of the handle to allow it to be squeezed in the conventional manner, or pushed using the palm of the hand.



Twin Gripping Jaws

The jaws have sculptured liners which give an excellent grip. The Pick-Up is operated by a wire not a cord eliminating fraying or stretching.

Magnetic Head



There is a magnet attached to the jaw which allows easy retrieval of pins, paper clips and other small ferrous items.

Rotatable Head/Stick Clip



The head can easily rotate to suit individual tasks. The Pick-Up is supplied with a removable stick clip, allowing it to be clipped on to wheelchairs and walking frames etc.

Part #	Description
PAT-AA8060Y	Standard Reacher
PAT-AA8063Y	Long Reacher
PAT-AA8064	Long with Lock and Forearm Support
PAT-AA8073	Extra Long



Handi-Reacher™

The Handi-Reacher™ is an ultra-lightweight and robust reacher that assists those with impaired reach and mobility. Featuring a single moving jaw operated with a trigger action from the handle, the Handi-Reacher™ comes with a magnet at the end for picking up small metal objects. The jaw rotates 360° to avoid the need to turn the wrist. The trigger — operated by four fingers instead of the typical two — is easy to grip and use.




Part #	Description	Length	Weight
PAT-AA8054Y	Standard	61 cm (24")	180 g (6.3 oz)
PAT-AA8056Y	Long	76 cm (30")	200 g (7 oz)
PAT-AA8058Y	Extra Long	90 cm (35½")	220 g (7.8 oz)



Combi Reacher

Lightweight and robust, this versatile new reacher is ideal for picking up everyday items around the home for those who have difficulty bending or reaching. Containing no skin traps, the comfortable four finger trigger action is suitable for those with a limited grip, whilst the 360° revolving head avoids the need to rotate the wrist, which may be painful for some users.

Part #	Description	Length
NEW-DAYRCHSTD	Standard	60 cm (24")
NEW-DAYRCHLONG	Long	80 cm (30")





Bed & Seating Equipment



Over Armchair Table

We've improved the design of our best selling Over Armchair tilt table. The Over Armchair Table has a curved, laminated teak tabletop that is easy to adjust in height and angle. The modern style includes a small table that can be fitted to the main table's left or right side. Drink holder fits left or right side. Improved fixation of tabletop.

- ARTG 230215

Specifications:

Table top	109 x 41 cm (43 x 16")
Side table top	20 x 20 cm (8 x 8")
Overall width	13.1 cm (5 1/2")
Height adjustment	75 - 91.5 cm (29 1/2 - 36")
Internal width	90.6 cm (35 1/2")
Maximum load	10 kg (22 lbs)



Part #	Description
HLC-758	Over Armchair Table



Tilting Overbed Table with Castors

Manufactured from chrome plated and powder coated steel with an adjustable laminated teak effect table top. The table is adjustable in height and comes with swivel castors for ease of manoeuvrability. There are rails on the edges of the tables to prevent items slipping off. Available in fixed or flat-packed

- ARTG 230215

Specifications:

Height adjustment	61 - 91 cm (24 - 36")
Table size	61 x 40 cm (24 x 15¾")
Castors	5 cm (2")
Maximum load	10 kg (22 lbs)



Flat-Packed
PAT-091558261



HLC-751C

Part #	Description
HLC-751C	Fixed
PAT-091558261	Flat-packed

Overbed Table

This is a great sturdy overbed adjustable height table, ideal for hospitals or in the home. Epoxy coated frame with a vinyl top. It also has recessed edges, which makes for easy cleaning. The 50mm twin wheel castors allows the table to be easily moved away from the bed when not in use.

Specifications:

Height adjustment	71 - 105 cm (28 - 43¼")
Tabletop	80 x 40 cm (27¾ x 14¾")
Height of caster wheels	5 cm (2")
Maximum load	15 kg (33 lbs)



Part #	Description
JAN-758KD	Overbed Table



Wooden Chair Raisers

Wooden Chair Raisers are cubes in bamboo material with three 1.6 cm (5/8") recesses of different diameters to suit various leg widths. Supplied in a set of four.

Specifications:

Height	9.3 cm (3¾")
Height raise	7.5 cm (3")
Maximum diameter of leg	3, 5 or 7.5 cm (1¼, 2 or 3")
Weight	425 g (15 oz) each
Maximum load	380 kg (60 st; 837½ lbs)

Part #	Description
PAT-AA3306	Wooden Chair Raisers - Set of 4



Stackable Cone Raisers — Bed or Chair

Stackable Cone Raisers are lightweight but strong moulded cones with a broad base for stability. Each features a recess for the chair or bed leg and anti-slip pads to hold the raiser in place. Set of four available for bed or chair.

Specifications:

	Chair	Bed
Height	9.3 cm (3¾")	14 cm (5½")
Height raise	7.5 cm (3")	7.5 cm (3")
Maximum diameter of leg	3, 5 or 7.5 cm (1¼, 2 or 3")	3, 5 or 7.5 cm (1¼, 2 or 3")
Weight	105 g (3.7 oz)	145 g (5.2 oz)
Maximum load	222 kg (35 st; 489½ lbs)	222 kg (35 st; 489½ lbs)

Part #	Description
PAT-AA3300Y	Stackable Cone Chair Raisers - Set of 4
PAT-AA3302	Stackable Cone Bed Raisers - Set of 4



Wooden Bed Raisers

Wooden Bed Raisers are black hardwood cubes in bamboo material with a recess in the top and a wooden base for more stability. Supplied in a set of four.

Specifications:

Height	17.8 cm (7")
Raise	13 cm (5")
Maximum leg diameter	7.5 cm (3")
Base	15.2 cm (6") square
Weight	750 g (1¾ lbs) each
Maximum load	380 kg (60 st; 837¾ lbs)

Part #	Description
PAT-AA3308	Wooden Bed Raisers - Set of 4



Adjustable Bed Cradle

Keeps covers off legs and feet. Can be used at the side or end of bed. Base slips between mattress and spring. Consists of a two-part frame. Flat for easy storage and traveling. Latex free

Specifications:

Height	45.7 - 68.6 cm (18 - 27")
Width	35.6 cm (14")
Base	71.1 cm (28")
Weight	3 kg (6.61 lbs)



Part #	Description
PAT-AA3652	Adjustable Bed Cradle

Bed Rail

The Bed Rail is designed to prevent the user from rolling out of bed. Once the user is in bed, simply slide the Bed Rail under the mattress. The weight of the user secures the Bed Rail. Constructed from sturdy epoxy coated steel..

Specifications:

Height	50 cm (19½")
Width (base)	58 cm (22¾")
Length (rail)	60 cm (23½")
Weight	2.5 kg (5.5 lbs)

*Now ordered in Box Qty

Part #	Description	Qty/Box
COR-AC38-4	Bed Rail	4



Bed Grab Rail with Legs

The Bed Grab Rail with Legs is ultra-sturdy, offering stability and confidence to those who have difficulty getting into and out of bed. The height-adjustable frame can be positioned to fit any divan or similar bed style. Two legs rest on the floor, and three colour-coded straps (supplied) secure the frame. Easy to assemble.

Specifications:

Handle height	85 - 95 cm (33½ - 37½") from the floor
Handle width at top	26.5 cm (10½")
Size of under-mattress frame	62 x 32.5 cm (24½ x 12¾")
Weight	3.5 kg (7¾ lbs)

127 kg
20 st 280 lbs
Maximum user weight

Part #	Description
PAT-AA3475	Bed Grab Rail with Legs



Rope Ladder Hoist

The Rope Ladder Hoist offers a simple way to assist with sitting up in bed. The bottom ends of the cords tie around the bed legs, and plastic ladder rungs enable the user to pull themselves up to a sitting position with a hand-over-hand movement. Suitable for single or double beds.

Specifications:

Rung length	20 cm (8")
Diameter	2.5 cm (1")
Cord length	3 m (118")
Weight	285 g (10 oz)
Maximum load	60 kg (9½ st; 132¼ lbs)

Part #	Description
PAT-AA3660	Rope Ladder Hoist



Adjustable Back Rest

The Adjustable Back Rest is secure and comfortable and will adjust to many positions, from low-lying angles to upright sitting positions. Breathable nylon fabric covers the metal frame. Includes a small pillow for additional comfort.

Specifications:

Maximum height of back	58 cm (22¾")
Width	65 cm (25½")
Depth	50 cm (19½")
Weight	2.8 kg (6¼ lbs)
Maximum load	113 kg (17¾ st; 249¼ lbs)

Part #	Description
PAT-AA3614	Adjustable Back Rest



Overbed Pole

This Overbed Pole can be used in conjunction with a conventional or hospital bed. It is lightweight, sturdy and easy to assemble. Also features a handle, which is on an adjustable webbing strap.

Specifications:

Height	168 cm (66")
Footprint	76 - 85 cm (30 - 33½")
Handle height	109 - 132 cm (43 - 52")
Weight	13 kg (28½ lbs)

135 kg
21¼ st 297½ lbs
Maximum user weight

Part #	Description
CZY-DAYOBPOLE	Overbed Pole



Plastic Curved Transfer Board

The Plastic Curved Transfer Board enables various positions to assist with functional transfers for individuals with limited strength and mobility. The transfer board is designed to aid seated transfers of individuals between two surfaces, such as bed to wheelchair or wheelchair to chair or toilet. The board helps reduce the risk of injury and strain on the caregiver assisting with transfers.

Specifications:

Length	71 cm (27½")
Width	37 cm (14½")
Weight	1.95 kg (4¼ lbs)



Part #	Description
PAT-AA8835	Plastic Curved Transfer Board

152 kg
23¾ st 335¼ lbs
Maximum user weight



Orbi-Turn Transfer Aid

The Orbi-Turn Transfer Aid features a steel frame and is designed for caregivers to facilitate standing and seat-to-seat transfers while minimising the risk of injury to the patient and themselves. The Orbi-Turn enables a patient with moderate upper and lower limb strength to participate in functional transfers while ensuring good working postures for caregivers. A height-adjustable, angled knee pad facilitates typical biomechanical hip and knee movement.

- ARTG 210114

Part #	Description	
PAT-091166883	Lowest position	Highest position
Bottom handle	81.5 cm (32")	117.5 cm (46¼")
Top handle	92 cm (36¼")	131.5 cm (51¾")
Knee pad (mid point)	23.5 cm (9¼ lbs)	51 cm (20")

190 kg
30 st 419 lbs
Maximum user weight



Standard Utility Chair

Adjustable height, contoured backrest and padded arm rests, vinyl upholstery and powder coated steel frame.

- AS/NZS 4688.2-2000
- ARTG 174748

Specifications:

External width	58 cm (22¾")
Seat depth	46 cm (18")
Seat width	45 cm (17¾")
Seat height	45 - 60.5 cm (17¾ - 23¾")
Back height	46.5 cm (18¼")



Part #	Description	Colour
ZHI-ZH026	Standard Utility Chair	Teal/Blue
ZHI-ZH026-MOCHA	Standard Utility Chair	Mocha



High Back Chair

High back chair, adjustable height, contoured backrest and padded arm rests, vinyl upholstery and powder coated steel frame.

- AS/NZS 4688.2-2000
- ARTG 174748

Specifications:

External width	64 cm (25")
Seat depth	48 cm (18¾")
Seat width	49 cm (19¼")
Seat height	45.5 - 60.5 cm (18 - 23¾")
Back height	73 cm (28¾")



Part #	Description	Colour
ZHI-ZH032	High Back Chair	Mocha
ZHI-ZH032-TEAL/BUE	High Back Chair	Teal/Blue
ZHI-ZH032CASTOR	Optional Rear Wheels - Set of 2	



Matching Legrest Stool

The adjustable height legrest stool is designed with a 10° angle for comfort and support and should be placed with the higher end against the chair. The legrest stool matches the high back chair (ZHI-ZH032) and is available in matching mocha or teal colours.

- AS/NZS 4688.2-2000

Specifications:

Seat pad	40 x 50 x 4 cm (15¾ x 19¾ x 1½")
Height	42-57 cm (16½-22½")



Part #	Description	Colour
ZHI-ZHSTOOL-TEAL	Legrest Stool	Teal
ZHI-ZHSTOOL-MOCHA	Legrest Stool	Mocha



Comfort Chair

The ultimate in comfort and mobility the Comfort Chair is an easy care solution for elderly or special needs seating. Cushioned adjustable air cell support has been incorporated into the design. The four lateral support wings fold away for easy transferring of the patient. This chair has a reclinable seat along with an adjustable backrest and calf rest.

Key features:

- Adjustable backrest (90°-170°) and calf rest
- Reclinable seat
- Air system for pressure relief
- Full stainless steel frame
- Three lockable castors and one steering castor (non-marking)
- Breathable polyurethane cover (fire retardant)
- Replaceable & easy to clean cover
- Four folding wings
- Standard dining tray
- Magazine pocket
- 10 years warranty on frame

Specifications:

Seat width	51 cm (20")
Seat depth	50 cm (19¾")
Seat height	63 cm (24¾")
Overall width	73 cm (28¾")
Overall length	135 - 188 cm (53 - 74")
Table to floor	88 cm (34½")
Table slide	23 cm (9")
Table size	40 x 40.5 x 2 cm (15¾ x 16 x ¾")
Back reclining	102°-136°
Seat tilting	0°-18°
Leg elevating	70° down to 5° up
Weight	43.5 kg (95¾ lbs)



Head Support Options



Part #	Description
ANI-DAYCOMCA1	Comfort Chair

Tilt Positions of Chair





Kitchen & Dining



E — complete set including Knife, Fork and Spoon

Queens Cutlery

Queens Cutlery includes a range of stainless-steel utensils with lightweight, ABS plastic built-up handles designed to help those with weak or restricted grips. Each handle is cylindrical with a slight taper away from the utensil and smooth contours.

- **Queens Standard Cutlery** is a popular range of standard utensils mounted in the Queens built-up handle. The knife, fork, spoon, junior spoon and teaspoon are available individually or as a set. Dishwasher-safe.
- **Queens Angled Cutlery** features stainless steel utensils with angled heads to help individuals with limited wrist or arm movement. The fork and spoon are available as left- or right-hand options. Dishwasher-safe.
- **One-handed Cutlery** features a stainless steel utensil with the Queens built-up handle and combines multiple cutlery functions. They are designed to help amputees and people with semi-paralysis after a stroke eat independently.

Specifications:

Handle length 9.8 cm (4⁵/₈"")
Handle diameter 3 cm (1¹/₄"")

Part #	Description	Weight
Standard Cutlery		
PAT-AA5500A	A — Knife	45 g (1.5 oz)
PAT-AA5501A	B — Fork	45 g (1.5 oz)
PAT-AA5502A	C — Spoon	45 g (1.5 oz)
PAT-AA5503A	D — Teaspoon	45 g (1.5 oz)
PAT-AA5517Y	E — Queens Retail Cutlery Set - One of Each	
Angled Cutlery		
PAT-AA5510A	Knife (not illustrated)	45 g (1.5 oz)
PAT-AA5511RA	F — Right-handed Fork	55 g (1.9 oz)
PAT-AA5511LA	G — Left-handed Fork	55 g (1.9 oz)
PAT-AA5512RA	H — Right-handed Spoon	55 g (1.9 oz)
PAT-AA5512LA	I — Left-handed Spoon	55 g (1.9 oz)
One-handed Cutlery		
PAT-AA5507A	J — Splayed Spork	45 g (1.5 oz)
PAT-AA5508A	K — Nelson Knife	55 g (1.9 oz)



Caring Cutlery

Caring Cutlery is a set of stainless steel utensils featuring contoured plastic handles with built-up ends. These utensils are ideal for individuals with weak grips and restricted wrist or finger movement. The knives and forks have a shaped indent on the top for the index finger to help with directional control. Also available as a four-piece set including knife, fork, spoon and teaspoon.

Specifications:

Handle length	12.7 cm (5")
Utensil weight	45-60 g (1.6-2.1 oz)

Part #	Description	Weight
Standard Cutlery		
PAT-AA5570	A — Knife, Ivory	50 g (1.8 oz)
PAT-AA5571	B — Fork, Ivory	50 g (1.8 oz)
PAT-AA5572	C — Spoon, Ivory	60 g (2.1 oz)
PAT-AA55721	D — Teaspoon, Ivory	45 g (1.6 oz)
PAT-AA55700Y	E — Cutlery Set - One of Each	215 g (7.6 oz)
Angled Cutlery		
PAT-AA5573R	F — Right-handed Spoon, Ivory	60 g (2.1 oz)
PAT-091098680	G — Left-handed Fork, Ivory	50 g (1.8 oz)
PAT-091098698	H — Right-handed Fork Ivory	50 g (1.8 oz)
Junior Cutlery		
PAT-91098706	Junior Knife, Sky Blue - Handle 102mm (4")	24 g (0.85 oz)
PAT-91098714	Junior Fork, Sky Blue - Handle 102mm (4")	22g (0.78 oz)
PAT-91098722	Junior Spoon, Sky Blue - Handle 102mm (4")	24g (0.85 oz)



E — complete set including Knife, Fork, Spoon and Teaspoon



Lightweight Foam-handled Cutlery

Lightweight Foam-handled Cutlery offers a comfortable design to help those with weak or restricted grips. The stainless steel utensils feature light grey foam handles that are ribbed and soft. Hand-washable only.

- **Standard Cutlery** including a knife, fork, spoon and teaspoon, is available individually or as a set. Hand-washable only.
- **Angled Cutlery** features stainless steel utensils with angled heads to help individuals with limited wrist or arm movement. The fork and spoon are at 40° angles and available as left- or right-hand options. The knife is at a 90° angle and can be used without cocking the wrist.
- **One-handed Cutlery** features a stainless steel utensil with the lightweight foam handle that combines multiple cutlery functions. They are designed to help amputees and people with semi-paralysis after a stroke eat independently.

Specifications:

Handle length	11.5 cm (4½")
Handle diameter	2.8 cm (1¼")



Part #	Description	Weight
Standard Cutlery		
PAT-AA5550	A — Knife	45 g (1.6 oz)
PAT-AA5551	B — Fork	34 g (1.2 oz)
PAT-AA5552	C — Spoon	42 g (1.5 oz)
PAT-AA5553	D — Teaspoon	30 g (1.1 oz)
PAT-AA5554	E — Foam Handle Cutlery Set - One of Each	151 g (5.3 oz)
Angled Cutlery		
PAT-AA5555	F — Knife	45 g (1.6 oz)
PAT-AA5556R	G — Right-handed Fork	40 g (1.4 oz)
PAT-AA5556L	H — Left-handed Fork	40 g (1.4 oz)
PAT-AA5557R	I — Right-handed Spoon	46 g (1.6 oz)
PAT-AA5557L	J — Left-handed Spoon	46 g (1.6 oz)



Newstead Cutlery

Newstead Cutlery features ergonomic, flexible, moulded plastic handles designed to adapt to the user's grip. With a comfortable rounded shape, it fits naturally into the palm and is ideal for people with arthritis, neurological impairments and other conditions that cause poor grip. All pieces, made of stainless steel and containing no latex, are dishwasher safe.

- **Standard cutlery set** includes a knife, fork, spoon and teaspoon.
- **Angled cutlery** features angled heads to help individuals with limited wrist or arm movement. The fork and spoon are at 40° angles, available as left- or right-hand options and may be bent to adjust the angle. The angled knife can be used for cutting without cocking the wrist.
- **Weighted cutlery** features standard cutlery heads but with a heavier handle, which helps reduce tremors to reduce spills.
- **One-handed cutlery** make dining with one hand easier for people with limited hand movement after a stroke, amputees and others who need to eat with one hand.

Specifications:

Standard		Weighted	
Weight	80 g (2.8 oz) per utensil	Weight	130 g (4.6 oz) per utensil
Set	320 g (11.3 oz)	Set	520 g (1¼ lbs)

Part #	Description	Weight
Standard Cutlery		
PAT-091097872	A — Newstead Cutlery Set	320 g (11.3 oz)
Angled Cutlery		
PAT-091097880	B — Right-handed Fork	80 g (2.8 oz)
PAT-091097898	C — Right-handed Spoon	80 g (2.8 oz)
PAT-091097906	D — Left-handed Fork	80 g (2.8 oz)
PAT-091097914	E — Left-handed Spoon	80 g (2.8 oz)
PAT-090197922	F — Angled Knife	80 g (2.8 oz)
Weighted Cutlery		
PAT-091097989	G — Weighted Cutlery Set	520 g (18.3 oz)
One-handed Cutlery		
PAT-091097997	H — Splayed Fork	80 g (2.8 oz)
PAT-091098003	I — Nelson Knife	80 g (2.8 oz)



A — complete set including Knife, Fork, Spoon and Teaspoon



G — complete weighted set including Knife, Fork, Spoon and Teaspoon



Sure Grip™ Cutlery

Comfortable ribbed built-up handles help people with weak grasp to hold utensils securely. Stainless steel with rubber handles. Dishwasher safe up to 52°C. Latex free.



Part #	Description
PAT-A703203	A — Knife
PAT-A703202	B — Fork
PAT-A703201	C — Table Spoon
PAT-A703200	D — Tea Spoon

Sure Grip™ Bendable Cutlery

The spoons and fork feature a metal shaft that is easily bent to any angle to accommodate reduced range of motion and to provide the greatest comfort for the user. (NOTE: the knife does not bend.) The rocker knife cuts with a minimum of motion. They have comfortable rubber built-up handles and stainless steel utensils. Dishwasher safe up to 52°C. Latex free.



Part #	Description	Length	Weight
PAT-A703208	A — Knife (does not bend)	210 mm	20 g (0.7 oz)
PAT-A703207	B — Fork	205 mm	21 g (0.74 oz)
PAT-A703206	C — Table Spoon	200 mm	21 g (0.74 oz)
PAT-A703205	D — Tea Spoon	82 mm	23 g (0.8 oz)

Homecraft Multi Holder

Comfortable, easily adjustable, strap handle that provides a secure grasp on a variety of handle shapes and sizes. Made from a soft, stretchy fabric, this utensil securely holds a range of items such as cutlery, pens, pencils, brushes and combs. Useful for people with limited hand function or a weak grip. One size fits all. For left or right handed use. Colour may vary. Not suitable for dishwashers.



Part #	Description	Length	Weight
PAT-AA5587	A — Splayed Fork	82 mm	22 g (0.78 oz)

Homecraft Nelson Knife

A stainless steel knife designed for onehanded eating. The stainless steel knife blade can be used with a rocker action and the pronged end used as a fork. A built up sleeve from the Kings Modular cutlery range may be added to the handle. Not suitable for industrial dishwasher temperatures. Handle length 100mm.



Part #	Description	Length	Weight
PAT-AA5597Y	Nelson Knife	100 mm	54 g (1.9 oz)




Kings Modular Cutlery

This attractive and comprehensive range allows a large amount of bespoke cutlery items to be made from 33 modular pieces. It caters for virtually all needs - economically, visually, aesthetically and functionally. The range covers: standard, angled and specialised utensils; built up, slim, contoured, lightweight and heavyweight handles; and a unique double ended feature to create special utensils for one handed eating or to provide an inconspicuous hand strap. All items are hygienic, non-corroding and dishwasher safe.

Kings Standard Built-Up Handle

The Standard Handle has a hollow cylindrical body. The Slim has an oval cross section and this handle is also used as the core for the built-up sleeves. The Contoured has a curved shape and finger contours. The Mushroom is a round hollow doorknob type handle which is easy to grip for those with weak or painful joints. The Heavyweight Handle is similar to the standard, but the core is filled with a metal bar to make it heavier. The Customised is custom-made by means of a metal inner core and pack of modelling material to make a unique personal shape.

Part #	Description	Width	Weight
PAT-AA5535	A — Contoured	22-30 mm	85 g (2.8 oz)
PAT-AA5536	B — Standard	25 mm	30 g (1 oz)
PAT-AA5539	C — Heavy Weight Handle	25 mm	156 g (5.5 oz)
PAT-AA5543	D — Slim	14-19 mm	29 g (1 oz)
PAT-AA5544	E — Mushroom	65 mm	52 g (1.8 oz)
Built-Up Sleeve			
PAT-AA5545	F — Oval Sleeve	23-30 mm	24 g (0.8 oz)
PAT-AA5546	G — Small Contoured	22 mm	19 g (0.7 oz)
PAT-AA5547	H — Large Contoured	29-35 mm	38 g (1.3 oz)



Kings Assessment Kit

This kit contains one of each utensil, handle and sleeve, plus the Kings Button Hook. Ideal for assessing individual patient needs. Supplied in an attractive case with a specific place for each item.



Part #	Description
PAT-AA5548	Complete Kit with Case

Kings Specialised Cutlery Utensils

The Splayed Fork (PAT-AA5527) is a combined fork and spoon and the Nelson Knife is a combined knife and fork. Both are designed to help with one-handed eating. The Small Rocker Knife is for cutting food with minimal rocking movement and can be used in the double-ended cutlery handles to make one utensil more versatile. The Large Rocker Knife has a broad semi-circular blade which is ideal for preparing larger food without cocking the wrist.



Part #	Description	Length	Weight
PAT-AA5527	A — Splayed Fork	82 mm	22 g (0.78 oz)
PAT-AA5528	B — Nelson Knife	110 mm	25 g (0.88 oz)
PAT-AA5529	C — Small Rocker	35 mm (width 60 mm)	19 g (0.7.5 oz)
PAT-AA5533	D — Large Rocker	121 mm - blade width	45 g (1.6 oz)

Kings Standard Cutlery Utensils

These standard utensils are mounted in plastic spigots that push tightly into the chosen handle but allow the blade to be angled as desired.



Part #	Description	Length	Weight
PAT-AA5520	A — Knife	112 mm	17 g (0.6 oz)
PAT-AA5521	B — Fork	104 mm	21 g (0.74 oz)
PAT-AA5522	C — Standard Spoon	97 mm	23 g (0.8 oz)
PAT-AA5523	D — Junior Spoon	91 mm	21 g (0.74 oz)
PAT-AA5524	E — Teaspoon	84 mm	17 g (0.6 oz)
PAT-AA5526	F — Soup Spoon	80 mm	21 g (0.74 oz)

Kings Angled Cutlery Utensils


Angled blades are ideal for people with limited wrist or arm movement. The knife is 52mm at 90° and can also be used as the second utensil in the double-ended handles. Forks and spoons have 40° angles. The angled adaptor adds extra length and a 30° angle to any utensil.



Part #	Description	Length	Weight
PAT-AA5530	A — Knife (8cm offset)	132 mm	20 g (0.7 oz)
PAT-AA5531R	B — Right Handed Fork	82 mm	21 g (0.74 oz)
PAT-AA5531L	C — Left Handed Fork	82 mm	21 g (0.74 oz)
PAT-AA5532R	D — Right Handed Spoon	82 mm	23 g (0.8 oz)
PAT-AA5532L	E — Left Handed Spoon	82 mm	23 g (0.8 oz)
PAT-AA5534	F — Adaptor	60 mm	10 g (0.35 oz)

Kings Lightweight Foam Handles


Made from grey sculptured closed cell foam for those unable to cope with the weight of standard built up handles. Washable.



Part #	Description	Width	Weight
PAT-AA5538	A — Small	20 mm	10 g (0.35 oz)
PAT-AA5537	B — Large	28 mm	12 g (0.42 oz)

Kings Soft-Coated Spoons

These spoons have been dipped in soft food quality PVC. This makes the spoons ideal for users with bite reflex problem when eating. Due to the type of use to which these item



Part #	Description	Length	Weight
PAT-AA5518	A — Spoon	100 mm	32 g (1.1 oz)
PAT-AA5519	B — Teaspoon	80 mm	25 g (0.88 oz)




Reflex Kitchen Knife Collection

Reflex Kitchen knives offer a secure but easy-to-grip handle, which is contoured and comfortable. The handle has been designed to counterbalance the weight of the knife. The angled, stainless steel blade minimises the user's wrist movement, making this collection ideal for those with weakened grip strength. Dishwasher-safe.

Specifications:

	Blade length	Weight
Reflex Slicing Knife	25.4 cm (10")	125 g (4.4 oz)
Reflex Carving Knife	20.3 cm (8")	155 g (5.5 oz)
Reflex Chef's Knife	20.3 cm (8")	165 g (5.8 oz)
Reflex Preparation Knife	12.7 cm (5")	140 g (4.9 oz)



A - Slicing Knife : A high quality stainless steel blade that is long and narrow. Ideal for slicing bread. Dishwasher safe.

B - Carving Knife : This high quality stainless steel blade with a scalloped rather than a smooth blade which will stay permanently sharp. Dishwasher safe.

C - Chef's Knife : This high quality stainless steel blade with a smooth sharp edge is useful for cutting or chopping food. Dishwasher safe.

D - Preparation Knife : This knife with a scalloped edge is useful for cutting and preparing food. The blade will not require sharpening. Dishwasher safe.


Part #	Description
PAT-AA5246Y	A — Slicing Knife
PAT-AA5244Y	B — Carving Knife
PAT-AA5245Y	C — Chef's Knife
PAT-AA5240Y	D — Preparation Knife

Reflex Comfort Grip Utensils

This range of Reflex Comfort Grip Kitchen Utensils has been designed with a soft-touch overmoulded handle to provide a more comfortable slip resistant grip. The handle is set at an angle to minimise wrist movement when using the product and the grip is slightly larger than the standard range to reduce the amount of finger movement required to grip the handle. The utensils are available with the same stainless steel blades as in the standard Reflex range. Dishwasher safe.

Specifications:

	Blade length	Weight
Reflex Carving Knife	20.5 cm (8")	214 g (7.5oz)
Reflex Chef's Knife	20.5 cm (8")	198 g (7 oz)
Reflex Preparation Knife	25.4 cm (10")	192 g (6.8 oz)
Reflex Slicing Knife	12.7 cm (5")	176 g (6.2 oz)



A - Chef's Knife : This high quality stainless steel blade with a smooth sharp edge which is useful for cutting or chopping food.

B - Carving Knife : This high quality stainless steel blade with a scalloped rather than a smooth blade which will stay permanently sharp.

C - Slicing Knife : A high quality stainless steel blade that is long and narrow. Ideal for slicing bread.

D - Preparation Knife : A high quality stainless steel blade that is long and narrow. Ideal for slicing bread.

Part #	Description
PAT-091207810	A — Chef's Knife
PAT-091207802	B — Carving Knife
PAT-091207828	C — Slicing Knife
PAT-091207794	D — Preparation Knife





Kitchen Workstation

The Kitchen Workstation is a multifunction food preparation unit that assists with many day-to-day kitchen tasks. The large, durable chopping board features a raised, L-shaped corner that can hold bread and toast in place while spreading, and removable stainless steel spikes secure bread, fruit or vegetables while being sliced, chopped or peeled. A removable clamp can hold various objects in place, such as a bowl, a tin or food, making food preparation easier. The removable grater and slicer mount to the chopping board, making one-handed use possible with minimum mess and waste. Four non-slip suction pads firmly secure the workstation in place. Dishwasher-safe.

Specifications:

Overall dimension	50 x 30 cm (19¾ x 11¾")
Clamp width	11 cm (4¼")
Weight	1.24 kg (2¾ lbs)



Part #	Description
PAT-AA5272	Kitchen Workstation
PAT-AA527201	Spare Clamp Set



Plate with Inside Edge

Textured bottom surface makes dish more attractive and durable. The special inside edge keeps food from sliding off the plate to make self-feeding easier for children and adults. Dishwasher safe to 80° and microwave safe. BPA, Phthalates and Latex free.

Specifications:


Diameter	22.86 cm (9")
Rim	2.4 cm (1")
Edge	1.27 cm (½")



Part #	Description	Colour
PAT-1425	Plate with Inside Edge	Off-White

Manoy Contoured Plates

These melamine plates are designed to help those with one good hand or who have difficulty picking up food with utensils. They are oval in shape with a sloping bottom and high sided end to help scoop up food without spilling it over the side. Not suitable for microwave ovens, but dishwasher safe up to 82°C.



Part #	Description	Size
PAT-AA5610	Manoy Contoured Plate Small	22.7 x 15.9 cm (9 x 6¼")
PAT-AA5612	Manoy Contoured Plate Large	27.9 x 19.7 cm (11 x 7¾")

Incurve Plate Surround

The Incurve Plate Surround is a flexible guard that can be fitted to plates to assist with one-handed eating. Three clips attach it to the rim and create a high inward sloping face to help with food collection and prevent spills. It is microwave- and dishwasher-safe.

Specifications:

External diameter	19 - 24.5 cm (7½ - 10")
Weight	45 g (1.6 oz)


Part #	Description
PAT-AA5662Y	Incurve Plate Surround

Hi-Lo Scoop Plate

Ideal for the one-handed user. The vertical wall guides food onto utensils. This lightweight plate has a non-slip rubbercoated base to prevent sliding. Entry wall steadily increases upward. This plate is microwave safe. Washes in top shelf of domestic dishwasher.

Specifications:

Diameter	23 cm (9")
Rim	1.9 cm (¾")
Wall	1.3 - 3.8 cm (½ - 1½")



Part #	Description	Colour
PAT-1540	Hi-Lo Scoop Plate	White

Round Scoop Dish

Moulded low in front and high in back, this unbreakable, round scoop dish is especially useful. Non-slip, rubber padded bottom provides control when using the scoop feature. Washes in top shelf of domestic dishwasher. Not recommended for microwave use. Latex free.

Specifications:

Diameter	20 cm (8¾")
----------	-------------



Part #	Description	Colour
PAT-081111780	Round Scoop Dish	Red
PAT-081005644	Round Scoop Dish	Ivory




Combination Cutting Board

A chef knife is attached so that it can be used in a chopping motion. Made of polypropylene. It is removable for cleaning and sharpening. The board is mounted on four suction cup feet. Three stainless steel spikes for securing foods are conveniently located on the board. Food guard aids with spreads. Dishwasher safe.

Specifications:

Base dimensions 40 x 30 cm (16 x 12")



Part #	Description
PAT-555662	Combination Cutting Board

Plastic Spreading and Chopping Board

The Spreading and Chopping Board is a durable, one-piece, moulded board with one raised, L-shaped corner, allowing jams and other foods to be spread on bread using one hand. It features six stainless steel spikes that can be inserted into the board and are designed to hold food while slicing or chopping.

Specifications:

Size 18.1 x 25.6 cm (7 1/8 x 10 1/8")
Weight 190 g (6.7 oz)



Part #	Description
PAT-AA528401	Plastic Spreading and Chopping Board

Bread Board

This hardwood board is fitted with nonslip feet, and six stainless steel spikes to grip bread or vegetables, when cutting one-handed. It also has an L-shaped lip in one corner to allow one-handed spreading of the cut slices. Size 325 x 175mm. Weight 972g.

Specifications:

Base dimensions 32.5 x 17.5 cm (12 3/4 x 6 3/4")




Part #	Description
PAT-AA5274	Bread Board

Plastic Spreading Board

The Plastic Spreading Board is a durable, one-piece moulded board with one raised, L-shaped corner, assisting the user when spreading jams and other foods on bread using one hand.

Specifications:

Size 18.1 x 25.6 cm (7 1/8 x 10 1/8")
Weight 190 g (6.7 oz)



Part #	Description
PAT-AA5284	Plastic Spreading Board



Caring Mug

The Caring Mug has a wide, stable base with two large contoured and angled handles, which enable it to be held in the correct position for drinking. The lid has a uniquely designed spout and matches the Caring Cutlery range. It is useful for those with tremors or who cannot sit upright. The mug is microwave- and dishwasher-safe. Supplied with a small aperture for liquids or an optional large aperture for semisolids.




Specifications:

Capacity 300 mL (10¼ fl oz)

Part #	Description
PAT-AA5710Y	Cup with Small Aperture

Shatterproof Mug

The 400-mL (13½ fl oz) Shatterproof Mug is made from a tough material to prevent accidental breakage. Liquid consumption is easily measured by marked gradations, and the mug is supplied with a choice of lids. The spout lid helps with drinking, and the recessed, spill-proof lid avoids splashes. The mug is microwave- and dishwasher-safe, but the lids should be hand-washed.



Part #	Description
PAT-AA5700	Shatterproof Mug

Handled Mug with Lids

The Handled Mug with Lids is a popular mug with large handles and shaped finger holds for ease of use. Liquid consumption is easily measured by marked gradations. The mug is dishwasher-safe, but the lids should be cleaned in warm, soapy water. Supplied with one spout lid and one splash lid.

Specifications:

Capacity 270 mL (9 fl oz)



Part #	Description
PAT-AA5720	Handled Mug with Lids

Feeding Cup

The Feeding Cup is designed for easy gripping and features a vented lid with a drinking spout. The cups are dishwasher-safe up to 135° C (275° F). Supplied in pairs.

Specifications:

Capacity 250 mL (9 fl oz)
Weight 80 g (3 oz)



Part #	Description
PAT-AA5750Y	4 mm Aperture - Set of 2
PAT-AA5752	8 mm Aperture - Set of 2

Insulated Mug with Lid

A polypropylene mug with an easy-to-grasp handle. Use for either hot or cold liquids. Lid regulates flow. Dishwasher safe to 82°C. BPA, Phthalates and Latex free. Not microwaveable. Capacity 355ml.



Part #	Description
PAT-1137	Insulated Mug with Lid

Cup with Temperature Regulated Lid

This version of a non-spill cup has a lid that changes colour when the contents are hotter than body temperature. If liquids hotter than 37.5°C (96°F) are put into the beaker the lid will change colour from blue to pink, alerting the user to be careful when drinking. The lid is safe and easy to use and can also be used with a straw. Dishwasher safe up to 130°C.



Part #	Description
PAT-AA5769	Cup with Temperature Regulated Lid

Delux Nosey Cup with Handles

This stylish cup is shaped to avoid contact with your nose when the cup is tipped allowing you to drink without tipping your head back. The material is tough polycarbonate which is shatterproof and resists cracking, prolonging the life of the cup. The ergonomic shape ensures that it sits comfortably in the hand, making it easy for those with limited grip. The cup is available with a set of handles which clip into the inset rim around the edge of the cup helping those with limited vision and allowing a firmer two handed grip if necessary. Simply align the small lugs on the cup to the corresponding holes in the handle and press firmly into place. To remove the handles, gently twist each side in turn and pull down. Microwave safe and dishwasher safe up to 80°C. Capacity 250ml.



Part #	Description
PAT-091536283	Delux Nosey Cup with Handles

Nosey Cutout Tumblers

These cups are redesigned for added comfort and convenience. A special cutout helps to maintain the proper head and neck positioning when swallowing, making drinking even easier. Each tumbler is transparent which enables at a glance monitoring of liquid intake. Durable polypropylene makes these tumblers dishwasher safe up to 109°C.



Part #	Description
PAT-1149	Nosey Cutout Tumbler 115 ml
PAT-1145	Nosey Cutout Tumbler 225 ml
PAT-1146	Nosey Cutout Tumbler 340 ml

Dysphagia Mug

Helps people who have trouble swallowing. Drinking while tucking the chin toward the chest gives normal swallowing mechanisms time to work. The design of this cup helps prevent liquids from escaping at the lips and directs the liquid to the centre of the mouth. It provides sufficient nose clearance to empty the cup without tilting the head back. The extended handle accommodates gripping by the entire hand or by the thumb only. Dishwasher safe to 82°C. Latex free.



Part #	Description
PAT-A779101	Dysphagia Mug Green 235 ml

One-Way Straws

These straws have one-way valves which stay filled with fluid even after removing the straw from a user's lips. This feature eliminates the possibility of sucking in too much air while drinking. Not suitable for dishwasher.



Part #	Description
PAT-081566314	One-Way Straws - Set of 10



Freehand Tray

The Freehand Tray is designed to enable one-handed carrying, leaving a hand to open doors, climb stairs and provide support. The unique handle folds flat inside the inner rim of the tray and lifts and locks into place in the upright position for carrying. The handle is easily removed for cleaning or for use as a standard tray. Dishwasher safe up to 130°C (266°F).

Specifications:

Size	46.5 cm x 30.5 cm (18¼ x 12")
Weight	1.1 kg (2½ lbs)
Max load	4 kg (8¾ lbs)



Part #	Description
PAT-AA5815	Freehand Tray
PAT-AA5810Y	Freehand Tray with Non Slip Mat

Stay Tray with Bean Bag

The Stay Tray with Bean Bag sits comfortably on your lap to provide a secure surface for eating, drinking or recreational activities, making it ideal for use on a chair or sitting in bed. The durable plastic tray has built-in handles and comes with a detachable bean bag base, which is removable and can be machine-washed. Includes a non-slip, heat-resistant silicon mat to stop plates or cups from sliding.

Specifications:

Size	46.5 cm x 30.5 cm (18¼ x 12")
Weight	1.1 kg (2½ lbs)
Max load	4 kg (8¾ lbs)



Part #	Description
PAT-AA5820	Stay Tray with Bean Bag

Breeze Indoor Rollator

Breeze Indoor Rollators offer a stylish alternative to a household trolley. Manufactured from painted aluminium tubing with four solid tyre wheels, each rollator includes a storage bag, locking brakes and height-adjustable handles.

Specifications:

Overall width	54 cm (21¼")
Overall depth	66 cm (26")
Handle height range	83.5 - 94.5 cm (34¾ - 37¼")



Part #	Description
PAT-091558360	Breeze Indoor Rollator

Adjustable Height Trolley Walker

The adjustable trolley walker is supplied with a cream painted steel frame. Clip on trays can easily be removed for cleaning. Adjustable height easily accessible loop lockable brakes. As the trolley has brakes, this Trolley can also be a walking aid.

- ISO 11199-2:2005
- ARTG 199708

Specifications:

Overall width	58.5 cm (23")
Overall depth	56 cm (22")
Handle height range	75.5 - 90.5 cm (29¾ - 35½")



*Now ordered in Box Qty

Part #	Description	Qty/Box
HLC-109-2	Beige Trolley Walker	2
HLC-109NS-2	White Trolley Walker	2
HLC-PL244	Replacement Tray	



Stainless Steel Cooking Basket

The Stainless Steel Cooking Basket enables users to cook vegetables without ever needing to lift a pan full of boiling water. Place vegetables into the basket, insert them into the pan, cook the food and remove the basket when complete. The large, flame-retardant, moulded nylon handle is easy to grip and offers heat insulation.

Specifications:

Top diameter	16.5 cm (6½")
Bottom diameter	12 cm (4½")
Depth	9 cm (3½")
Weight	245 g (8.6 oz)



Part #	Description
PAT-AA5351Y	Stainless Steel Cooking Basket



Homecraft Undo-It Jar and Bottle Opener

A V-shaped coated metal opener with non-slip lining that is screwed underneath a shelf or cupboard. The jar is slid into the V until the lid wedges. It is then twisted to remove the lid. Suitable for tops 25 to 75mm diameter. Weight 145g. 140 x 100mm. Fixings not included.

Part #	Description
PAT-AA5010Y	Undo-It Jar and Bottle Opener



Homecraft Twister Jar Opener

A cone shaped rubber moulding with fluted finger grips on the outside, and ribbing on the inside to give a good grip for opening twist-off jar lids. Suitable for tops 20 to 85mm diameter. Weight 100g.

Part #	Description
PAT-AA5000Y	Twister Jar Opener



Everyday Bib

The Everyday Bib is practical for everyday use. The colourfast polycotton outside is durable but absorbent. The soft and flexible waterproof backing helps reduce clothes stains. All bibs are fitted with press studs, so the user can form a pocket to catch excess food. Users can fasten the bibs with hook-and-loop patches for quick, easy donning or metal poppers for extra security, with two alternative fastening positions for the poppers. The bibs are fully washable and can be tumble-dried at low temperatures. Available in three sizes and two colours.

Specifications:

Small	58 cm x 38 cm (23 x 15")
Medium	68 cm x 46 cm (26¾ x 18")
Large	89 cm x 46 cm (35 x 18")

Part #	Description
PAT-091421148	Everyday Bib, Red, Small
PAT-091421155	Everyday Bib, Blue, Small
PAT-091099399	Everyday Bib, Blue, Medium
PAT-091099407	Everyday Bib, Red, Medium
PAT-091099365	Everyday Bib, Blue, Large
PAT-091099373	Everyday Bib, Red, Large





Homecraft is a registered trademark of Performance Health and/or its subsidiaries and may be registered in the United States and other countries. © 2022 All rights reserved. Unauthorised use is strictly prohibited.
HC22CAT



1300 473 422 | sales@performhealth.com.au | performhealth.com.au

Ausmedic Australia Pty. Ltd. trading as Performance Health ANZ
Unit 3, 3 Basalt Road Pemulwuy NSW 2145 Australia