

cura¹
Fall Prevention Products
Home Care



Our LED Pager is standalone and works with up to 3 devices making it ideal for home care.

Easy to use and economical.

3508 LED Pager

Features:




- *Can be programmed to monitor up to 3 wireless devices
- * Vibration and tone selection for hard of hearing.
- *Auto-reset when paired with a 2719 transmitter .
- *Alarm memory and auto re-arm
- *Direct signalling from wireless devices for use in home care environments.
- *Easily cleanable buttons for a hospital setting.
- *Power requirements: 2 x AA batteries (included).
- *Housing material: Durable ABS plastic

Performance Health
1300 4 REHAB (1300 4 73422)
sales@performhealth.com.au

HOW TO ORDER:

Start with the 3508 LED Pager and then choose up to 3 wireless programmable transmitters from the range listed below.

Picture	Code	Description
	3508	LED Pager



Wireless Programmable Transmitters: Choose the devices to best suit requirements. You can choose up to 3 devices. Add devices at any stage. All pads require a 2719 or 2720 transmitter.

Picture	Code	Description	Best Feature	Transmitter Required
	2630	Wireless Chair Pad	Chair pad can be placed anywhere in the room	2719
	2691	Universal Chair Pad	18mth battery life/warranty	2720
	2629	Wireless Above Mattress Bed Pad	Can detect people 40kg & under	2719
	3515	Wireless Under Mattress Bed Pad	People won't even know it is there	2719
	2670	Wireless PreAlert Bed Pad	Alerts when they move to the edge of the bed	-
	2690	Universal Above Mattress Bed Pad	18mth battery life/warranty	2720
	2624	Wireless Floor Sensor Mat	Can be moved around the room as required	2719
	2671	Universal Single Fold Floor Sensor Mat	18mth battery life/warranty	2720
	2672	Universal Dual Fold Floor Sensor Mat	18mth battery life/warranty	2720
	2692	Universal Floor Sensor Mat	18mth battery life/warranty	2720
	2719	Replacement Transmitter	12-15mth battery life	-
	2720	Replacement Universal Transmitter	18mth battery life/warranty	-
	2721	PIR Motion Sensor	Ideal as an out of/into room sensor	-
	3503	Door Reed Switch	Ideal as an out of/into room sensor	-
	2588	Pendant with Lanyard	Economical & lightweight	-
	2595	Pendant on Wristband	Low battery warning LED	-
	2599	SOS Button	Large button, easy to see and use	-
	3502	Handset with Reset	Easy to use wireless handset	-
	2558	Waterproof Call Point with Reset	Pull string & push button activation	-

cura1 GURANTEE: All Products supplied with full warranty and guaranteed to be free from defects in material and workmanship. The range of wireless devices may vary and be dependent on your specific environmental conditions. cura1 products are designed to augment well designed care programs and practices, cura1 does not claim that any of its products are suitable as a substitute, or as a total solution for individuals that require care and constant supervision, and are at risk.