

# Three Channel Cordless Falls Monitor



PART NUMBER: 2702

## FEATURES

- Modern ergonomic design
- Monitor up to 3 cordless devices simultaneously
- NO / NC switch in battery compartment
- Comes with batteries
- Alarm auto re-arm
- Amber LED visual alarm indicator
- Battery operated or optional external plug pack (2406)
- Up to 6 months battery life (nominal)
- Can be connected to a nurse call system
- 3 selectable volume levels
- Alerts for lost nurse call connection, lost Pad, lost AC power supply & low battery
- Alerts for low battery in paired bed / chair Pad transmitters

## TECHNICAL SPECIFICATIONS

- DC operation using either a 3 x AA batteries (included) or 6 V plug pack (2406)
- Dimensions 120 x 83 x 35 mm
- Weight 210 g without batteries
- Up to 50 m\* range when used with cordless bed pads and chair pads
- Operating temperature 5-40° C, relative humidity < 85%
- Operating frequency 433 Mhz fully supervised
- Volt free normally open relay contact alarm output
- Low battery detection
- Missing AC / DC power supply detection and indication
- Housing: ABS (UL-94HB)

## ACCESSORIES

**All Cordless Bed & Chair Pads & Floor Mats, 2721 iBeam, 3503 Wireless Door Contact, 2595 Wristband Alarm Transmitter, 2599 Wireless SOS Call Button, 3502 Wireless Handset, 2558 Wireless Call button with Reset & Pull Cord, 2588 PET, 2742 Soft Monitor Cover, 3504 ActiveCare Mounting Bracket, 2643 Nurse Call Cable, 2406 Plug Pack 6 Vdc**

cura<sup>1</sup>.com

\*The range of wireless devices may vary and be dependent on your specific environmental conditions.  
cura<sup>1</sup> products are designed to augment well designed care programs and practices, cura<sup>1</sup> does not claim that any of its products are suitable as a substitute, or as a total solution for individuals that require care and constant supervision, and are at risk.