

Standard Falls Monitor with Safetime



PART NUMBER: 2701

FEATURES

- Modern ergonomic design
- Pause mode to attend to a patient without triggering alarm
- Tamper-proof reset option using magnetic caregiver key
- Choice of tone, and volume control for alarm
- 0, 2 sec. and 15 min. alarm delay settings
- Port to connect with nurse call system (NC / NO)
- Operates on batteries or AC adaptor
- Low battery notification
- Input lost notification (Pad Lost, Nurse Call, AC lost)
- Change Pad (SafeTime) indication
- Comes with mounting bracket, hanger strap, soft cover, reset key and batteries

TECHNICAL SPECIFICATIONS

- Power supply: 6 Vdc, 300 mA & / or 3 x AA alkaline batteries.
- Battery life: 1 year (nominal)
- Operating temp: 5°C to 40°C
- Rel. humidity: < 85%
- Housing material: ABS (UL-94HB, flammability grade)
- Dimensions: 120 x 77 x 34 mm
- Net weight: 140 g (without batteries), 200 g (with cover)
- Pad / Mat socket: RJ10
- Adjustable volume: up to 108 db @ 10 cm
- Incontinence (Wetness): Under development

ACCESSORIES

- All Premium, Universal and SafeTime Bed & Chair Pads & Floor Mats.**
2406 Plug Pack 6 Vdc, **2642** Caregiver Reset Key,
3504 ActiveCare Mounting Bracket, **2636** Wheelchair / Chair Seat Belt, **2643** Nurse Call Cable.

cura¹.com